




PAL RSI and PAL RTC Sample Prep and Injection





Content

Introduction	4 - 7
PAL RTC - Robotic Tool Change	8 - 9
PAL RSI - Robotic Sample Injection	10 - 11
PAL DHR Dual Head	12 - 13
Tools and Modules	14 - 23
PAL Sample Control	24 - 25
PAL Method Composer	26 - 27
Applications	28 - 33
PAL Accessories & Consumables	34 - 43
Specifications	44 - 45



PAL System is the most widely used sample preparation and handling platform.

PAL System is your tool box for sample preparation, from a simple liquid injection to complete workflows. A PAL System can be adapted or extended to meet almost any requirement.

Numerous options allow to increase sample capacity or add further modules for sample preparation. Listed below are the method steps that a PAL System can perform:


- Liquid injection GC and LC
- Gas injection
- Headspace sampling
- Dynamic headspace sampling
- SPME and SPME Arrow sampling
- Temperature controlled storage of samples
- Incubation 35 200 °C
- Tool change (PAL RTC only)
- Transport of vials and other objects
- Vortex mixing
- Centrifugation up to 5000 g
- Dilution
- Standard addition
- Liquid / liquid extraction (LLE)
- Derivatizations
- μ Solid phase extraction (μSPE)

More than 50'000 users in gas and liquid chromatography, mass spectrometry and optical spectroscopy can't be wrong. Read about their success stories in "Ingenious News" where you get regular application updates from users.



Sample preparation is the key to success to achieve precise and accurate results.

For a typical chromatographical analysis 61 % of time is spent on sample prep, 30 % of errors are linked to sample prep*.



PAL RTC - Highest Productivity & Flexibility

The PAL RTC with Robotic Tool Change is the logical (r)evolution of the successful PAL-xt product line. It is a robotic platform for the efficient and safe automation of most sample preparation steps.

The robotic change of tools enables unattended 24/7 operation, even for multistep workflows and thereby greatly increases the productivity of labs. At the same time process safety is optimized since all operations become traceable.

PAL RSI - The Workhorse

The RSI is the robust workhorse for analytical labs. It is the cost effective solution for labs in need of robustness and uncompromising performance.

The PAL RSI is a safe investment: if your requirements grow an upgrade to the full RTC functionality is possible and gives access to complete workflow automation.

PAL Sample Control / Chronos Software for Efficient Operation of PAL RTC & RSI



The latest generation of PALs can all be controlled by the user friendly PAL Sample Control / Chronos software. It interfaces seamlessly with many common CDS and MS-data systems (e.g. Agilent Chemstation, Masshunter, Sciex Analyst, or Thermo Scientific Xcalibur) or can be used for offline sample preparation. With a few clicks you can import or generate sample lists and start the data acquisition. Or you can quickly set up workflows to eliminate tedious manual operations. PAL Sample Control / Chronos allows overlapping of time consuming steps. This increases sample throughput greatly and boosts productivity.

PAL Method Composer



For PAL Systems integrated in a GC-MS or LC-MS system, PAL Method Composer lets you easily create new methods.

The graphical user interface allows the creation of a method by drag & drop. The check for validity of the method is performed automatically on the fly. More details on page 27 or on PAL Method Composer



The PAL RTC is all about increasing productivity.

Robotic Tool Change takes productivity to a new level.





Automation is the way to increase productivity and (process) safety in the laboratory. Transferring repetitive or dangerous manual tasks to a robot improves safety. The possibility to run the instrument 24h/day increases throughput, especially for long sample preparations.

The PAL RTC was developed to maximize productivity in analytical and clinical labs. Robotic Tool Change (RTC) brings sample preparation to a higher level.

Every process requires a number of different tools for best performance, e.g. a 10 μ L syringe for the accurate addition of small volumes followed by the dilution with a 1mL syringe. Robotic Tool Change allows to switch between different tools automatically.

This additional versatility in combination with the large number of available tools enable the design of tailored automation processes.



PAL RTC with standard x-axis for GC & GC/MS

Ingeniously productive.

- Automatic selection of the syringe with optimal accuracy for adding standards or preparing serial dilutions
- On the fly switching between a syringe tool for the addition of an internal standard and the LC/MS Tool for subsequent analysis
- Possibility to permanently configure several workflows on one system for a walk-up prep station, e.g. Liquid/ Liquid Extraction and Solid Phase Extraction (SPE)
- Automated optimization of methods e.g. by selecting the most suitable Solid Phase Micro Extraction (SPME) fiber from an array of 4 different ones
- Derivatization reactions performed without manual intervention for productivity, protection against hazardous chemicals and process safety
- Automation of labor intensive manual workflows like protein digestion

For detailed examples of workflows see p. 28 - 33.

The PAL RSI is the reliable workhorse for analytical labs.

The PAL RSI defines the industry standard for intelligent sample preparation.



PAL RSI with compact x-axis for LC/MS

Ingeniously reliable.

Laboratories work under time pressure and often with a tremendous workload. The reliability of hardware and software should not be something the user has to worry about. Reliability is just expected from every tool.

That is exactly what the PAL RSI was designed for. It is a tool that you can rely on. 50.000 PAL systems worldwide are proof of this.




PAL RSI with standard x-axis for GC & GC/MS

The most flexible system on the market.

RSI is a flexible tool. Its open and modular architecture makes it the most versatile system on the market. Tools can be exchanged readily within minutes.

PAL customers working with GC love the possibility to use liquid, headspace or SPME sampling on the same system. LC customers use the PAL because of its huge sample capacity, the range of syringes and valves available.




PAL DHR Dual Head Productivity x2

- Two heads can move independently
- The heads can execute two independent workflows
- Headspace sampling with standard addition
- Parallel injections into two GC- or LC-systems
- The software optimizes the time table to maximize productivity
- Collision control → process safety and easy programming
- 80 to 200 cm wide working area
- Combinations of RTC/RSI heads
- With an RTC head toolchange expands the options further

Application Examples

- Automated weighing of samples, combined with the analysis of a weighed sample
- High troughput LCMS injections with 2 completely independent LC streams (2 injection valves, 2 columns or staggered injections into mutltiple streams.)
- Sample prep combined with parallel headspace GC analysis



Get weighed Injection using Data acquisition

For precise sample preparation you need the right tools and modules.
Here is the complete toolbox.



The right tool for each application

Starting from 10 mL of a water sample, extracting pesticides by SPE and precisely injecting 1 µL of extract: the PAL Sytem performs this job flawlessly. This is only one example how PAL System improves productivity and process safety.

Liquid, gaseous or solid samples require different workflows and therefore different tools. The following pages introduce you to the PAL Toolbox equipped with the best tools and modules available for sample preparation.

Change the tool at any time with the Park Station (PAL RTC only)

The unique Park Station allows a robotic tool change (syringes with different volumes or different tool types) for advanced sample preparation, liquid handling (dilutions), derivatization steps or any other time consuming repetitive step. The additional versatility in combination with the increased volume range are significant benefits and allow the realization of tailor-made automation processes.

- Park Station for up to three injection tools
- Installation of up to 2 Park Stations on one PAL RTC
- Process safety through tool recognition and position control
- Designed for automated screening of method conditions with multiple tools



PAL Tools

Liquid Syringe Tool

- Dedicated Tools for 57 mm or 85 mm syringe needle lengths available
- Syringes available:
- 1.2 μ L / 5 μ L / 10 μ L / 100 μ L with D7 tool
- 250 μL / 500 μL / 1'000 μL with D8 tool
- 5'000 μL / 10'000 μL with D18 tool



µSPE Tool

- Liquid Tool for handling μSPE cartridges
- Holds syringes from 100-1000 μL



SPME Fiber Tool

- SPME holder for increased fiber protection
- Compatible with a variety of SPME fibers
- Easy fiber exchange by hand
- For 10 mm or 20 mm fiber length
- SPME Fibers are available with the following sorbents:
- PDMS, DVB, PA, Carbon WR, PDMS/Carbon WR, Carbon WR/DVB



SPMF Arrow Tool

- For SPME Arrow fibers with 1.1 and 1.5 mm diameter
- Easy fiber exchange by hand
- SPME Arrows are available with the following sorbents:
- PDMS, DVB, PA, Carbon WR, PDMS/Carbon WR, PDMS/Carbon WR/DVB



Headspace Tool

- Dedicated tools for different syringe types available:
- 1 mL / 2.5 mL / 5 mL (with corresponding tools)
- Syringe temperature 40 °C up to 150 °C in 1 °C steps
- Syringe flush with inert gas flow through X-Y-Z rail
- Magnetic vial transport for 2 mL, 10 mL and 20 mL vials



ITEX DHS Tool

- For powerful dynamic headspace extractions with Tenax TA trap, other materials available
- Active cooling to minimize idle times
- Temperature range: 40 350 °C for efficient desorption



LCMS Tool

- The special design minimizes carryover even for the most sensitive detectors
- Enables special injection techniques like sandwiching a sample between air gaps for accurate small volume injections
- Flow control guarantees bubble-free solvent delivery
- Position detection for zero dead volume positioning of syringe needle in the injector port
- Injection volume LCMS-100, 1-80 μL
- Injection volume LCMS-250, 1-230 μL
- Needle length 57 or 80 mm





Pipette Tool

- For automated pipetting with 20 μL, 200 μL or 1000 μL tips
- Single or multi-dispense mode
- A special adapter allows direct injection into LC-valves
- Also works with capped vials in combination with the DeCapper Module





Gripper Tool

- For transporting non-magnetic items
- Different gripping adapters available
- 4-pin gripping adapter included for items 5-24 mm in diameter



Dilutor Tool

- For the addition of larger amounts of liquids
- With the special "transfer" mode also small volumes of liquids can be handled with high precision
- the Dilutor Tool with 80 mm needle allows aspiration of samples from 10/20 mL vials



PAL Modules

Agitator Module

- For the incubation and agitation of samples
- 6 positions for 20 mL vials
- Temperature range 40 200 °C
- Agitation speed 250-750 rpm
- Optional adapters for 2 mL or 10 mL vials



Barcode Reader Module

This unique Barcode Reader allows PAL RSI and PAL RTC to read the barcode labels on 2 mL, 10 mL and 20 mL vials regardless of the orientation on the vials. Therefore it ensures highest process safety and traceability.

- Reads horizontal 1D barcodes
- Two scanners allow identification of vials irrespective of
- position of the barcode
- Works with 2 mL / 10 mL / 20 mL vials



Centrifuge Module

- For the efficient centrifugation in automated workflows (e.g. protein precipiation, phase separation)
- Full software control with PAL Sample Control
- Control through many Chromatographic Data Systems
- Relative centrifugal force:
- 2000 x g with Centrifuge Combi
- 5600 x g with Centrifuge 2 mL
- 2600 x g with Centrifuge 10 mL
- Vial format:

Centrifuge Combi: 2 mL (4 vials, 10 mL / 20 mL (2 vials)

Centrifuge 2 mL: 2 mL (8 vials)

Centrifuge 10 mL: 10 mL (6 vials, optional 2 mL and 4 mL)









Rotor 2 mL

Rotor10 mL

Rotor Combi

19

DeCapper Module

- Opens/closes 2, 10, 20 mL screw cap vials without any change of hardware (no adapters required)
- Defined torque guarantees the reproducible and leak-tight closing of headspace vials



Dilutor Module

- For the efficient and accurate addition of larger amounts of liquids
- 100 μL, 1 mL, 5 mL & 10 mL dispensing syringes available
- Optional selector valves allow dispensing of up to 5 different liquids





Heatex Stirrer Module

- For powerful mixing and heating in sample prep and SPME Arrow
- Temperature range 30 ° 150 °C
- Stirring speed up to 1600 rpm (i.e. 200 cycloidal loops)
- Optimized for 20 mL vials (for 10 mL vials special adapters are required)



SPME Arrow Conditioning Module

- For the conditioning of SPME Arrow and SPME fibers prior to sample enrichment, max. 350 °C
- Position for automated conditioning
- Position for manual pre-conditioning
- Automated purge gas valve
- Manual gas valve for pre-conditioning



Vortex Mixer Module

- For efficient mixing (dilution / extraction).
- Standard vial sizes: 2 mL / 10 mL / 20 mL
- 1 additional slot for custom specific vials
- Provides efficient mixing with up to 2000 rpm



Standard Wash Module

- Wash Module for low volume injections
- 4x 10 mL wash solvent vials
- 1x 10 mL waste vial
- Optional Waste Port Adapter to connect a tube to a waste bottle



Fast Wash Module

- Cleans syringes of gauges 19 to 26
- Integrated pumps for active wash solvent delivery
- Supports two different wash solvents (aqueous and organic)
- Reduces wash solvent consumption by automatic flow adjustments
- Can be mounted underneath Valve Drives to minimize required space

Fast Wash HF Module

Wash Module for flow rates up to 40 mL/min

LCMS Wash Module

Passive Wash Module for use with LCMS-Tools or Dilutor Tools

Large Wash Module

- Wash Module for large volume injections
- 2 x 100 mL solvent container (glass) with septum cap
- Waste port with tubing olive to connect waste bottle



Solvent Module

- For large solvent demands, e.g. for the addition of larger volumes of solvents or serial dilutions
- 3 x 100 mL solvent container (glass) with septum cap



Tray Cooler Module

For the storage of 3 racks or plates under defined temperature conditions between 4 °C and 40 °C.

Capacity:

- 3x MTP (Multi Titer Plate)
- or 3x DW (Deep Well Plate)
- or 3x VT15 (15 x 10 mL)
- or 3x VT54 (54 x 2 mL)
- or 3x VT70 (70 x 1 mL)
- or combinations
- Depending on tray type inserts are required



LigCooler Module

- Trayholder incl. sample tray
- for up to 32 vials of 10/20 mL
- External circulator bath for heating/cooling not included.



Trayholder Module

The Trayholder offers sample storage at room temperature. A PAL with extended x-axis length can hold up to 7 Trayholders.

Capacity:

- 3x MTP (Multi Titer Plate)
- or 3x DW (Deep Well Plate)
- or 3x VT15 (15 x 10/20 mL
- or 3x VT54 (54 x 2 mL)
- or 3x VT70 (70 x 1 mL)
- or combinations
- or 60x 10/20 mL (with one R60 tray)



Flow Cell Module

- Flow cell for sampling from a flow-through stream
- Up to 6 flow cells on one holder



Stack Modules (6 DW and 12 MT)

For the storage of racks or plates at room temperature. Allows to store samples with light sensitive compounds. A maximum of four Stacks can be configured on a PAL with extended x-axis length.

Capacity 6 DW: Capaci

- 6x MTP (Multi Titer Plate)
- or 6x DW (Deep Well Plate)
- or 6x VT15 (15 x 10 mL)
- or 6x VT54 (54 x 2 mL)or 6x VT70 (70 x 1 mL)
- or combinations

Capacity 12 MT:

 12x MTP (Multi Titer Plate, shallow well)

Peltier Stack Modules (2DW, 6DW and 12MT)

For the storage of two (2DW) or 6 (6DW) racks or plates or 12MT plates (12 MT) under defined temperature conditions between 4 °C and 40 °C. Allows to use transparent standard vials with light sensitive compounds.

Capacity 2DW:

- 2x MTP (Multi Titer Plate)
- or 2x DW (Deep Well Plate)
- or 2x VT15 (15 x 10 mL)
- or 2x VT 54 (54 x 2 mL)
- or 2x VT70 (70 x 1 mL)
- or combinations

Capacity 12MT:

- Capacity 6DW:
- 6x MTP (Multi Titer Plate)
- or 6x DW (Deep Well Plate)
- or 6x VT15 (15 x 10 mL)
- or 6x VT 54 (54 x 2 mL)
- or 6x VT70 (70 x 1 mL)

or combinations

- 12x MTP (Multi Titer Plate, shallow well)

(Peltier Stack 6DW shown)

(Stack 6DW shown)

Valve Drive Module

Universal Valve Drive for applications like sample injection, column switching for online LC-LC/MS or online SPE-LC/MS, Multiplexing, column selection and many more.

- Valve Drive supports VICI/Valco and Rheodyne valve types
- Injection port bottom sensing minimizes carryover
- Constant Force Technology to reduce dead volume during
- Injection process
- UHPLC/HPLC: up to 50% faster switching times for optimum
- System performance and prolonged column life time
- Stackable design to reduce the space required
- Huge flexibility to arrange multiple valve solutions
- The Fast Wash Module is also stackable below a Valve Drive



Two Injection Valves Staggered configuration reduces space required.


MHE Module

• Multiple Headspace Extraction Module



For more information about PAL Tools and Modules please visit our website: www.palsystem.com





Task Description 1 Task Description 2 Wather-fapord Plat sample in Agitator 3 bight-or-fapord Repeated Repea

Powerful Method Editor

A powerful yet easy to use tool to generate custom methods is included.

PAL Sample Control / Chronos

CHRONOS



Easy to use routine software

PAL Sample Control / Chronos software is the user-friendly tool for the daily routine jobs. With a few clicks sample lists are generated or imported. Now PAL Sample Control / Chronos starts the operation and the data acquisition. Since PAL Sample Control / Chronos interfaces seamlessly with most of the major chromatographic or MS data systems only one sample list has to be handled. Different user levels ensure process safety.

Productivity

PAL Sample Control / Chronos allows overlapping of time consuming steps. It optimizes automatically the timing of various steps in a sample preparation process and generates a schedule that minimizes the runtimes of sequences. This increases sample throughput greatly and boosts productivity.

Powerful Method Editor

While PAL Sample Control / Chronos is straightforward to use in the daily operations it is also a powerful tool for the generation of tailored methods. A set of tested methods that comes with every system (e.g. headspace injection, partial loop liquid injection) can be used as templates and optimized or tailored for specific workflows. Furthermore a large number of building blocks (tasks) for method development are part of the software. These building blocks make it easy to generate new methods, even for complex workflows. PAL Sample Control / Chronos is required for operating the PAL DHR.

Supported CDS and MS-Data systems

- Agilent ChemStation, GC, LC and MSD
- Agilent MassHunter B07, B06 (GC)
- Agilent MassHunter LC
- Bruker MS Workstation 8.2
- Chromperfect 6.0.1 and higher
- DataApex Clarity 3.0 and higher
- EZChrom 3.21 and higher
- GL Sciences Evolution Workstation
- OpenLab ChemStation C01.06 and C01.07
- Sciex Analyst 1.41 and higher
- Shimadzu LabSolutions
- Thermo Xcalibur 1.4 bis 4.1.5.0, LC und GC
- Thermo Omnic 9
- Waters MassLynx 4.1




PAL Method Composer, methods in minutes




PAL Method Composer is a tool to intuitively generate methods for the integration in Chromatographic Data Systems (CDS). By simply dragging and dropping of the individual prep steps you can build a method in minutes. Each step's functions and parameters are explained. The instrument configuration is visualized to easily optimize the sequence.

PAL Method Composer gives users easy access to tailored methods.



The parameters of the steps are default values that were experimentally determined. However, each step can be adjusted for specific methods.



PAL Method Composer works with the following CDS:

- Agilent: ChemStation / MassHunter / OpenLAB
- Bruker: Compass HyStar
- Sciex: Analyst AAO and ADD
- Shimadzu: GCMSsolution Software / LabSolutions
- Thermo: Xcalibur / Chromeleon / TraceFinder
- Waters: Empower 3

Automated sample prep is reliable, traceable, productive.

The workflows and application examples (following pages) were all realized with PAL Sample Control.

The tools and modules applied are shown.

Processing of blood samples with protein precipitatation, directly from primary tubes

Aspirate 50 µL whole blood from

Dispense blood into 2 mL vial

Clean tool

Add 250 µL water

Add 750 µL precipitation solvent

Vortex 90 sec

Centrifuge vial @ 3000 g for 300 sec

Inject 5 µL of supernatant

Start LC-MS

Clean tool

Next sample

Liquid-liquid extraction (LLE), often used, works best if automated

Aspirate sample

Dispense into 10 mL vial

Clean tool

Add 2 mL water

Add 2 mL hexane

Vortex 60 sec

Inject 5 µL of supernatant

Clean tool

Next sample

Automated analysis of fatty acid methyl esters (FAME)

Start Trans-esterification

Add ISTD

Add NaOMe

Vortex vial

Wait (Derivatization time)

Add heptane

Vortex vial

Add citrate

Vortex vial


Inject sample

Start analytical run

End derivatization


Start Trans-esterification of next sample

The trans-esterification of fatty acids to generate the corresponding methyl esters (FAME) is a very common and at the same time tedious procedure. Its automation increases productivity and prevents exposure of humans to hazardous chemicals.



μ Solid Phase Extraction (SPE), the express way to clean samples

A wide variety of materials allows the selective clean-up of samples for further analysis.



You find all our application notes on:

Application Notes




Contact the experts for sample prep:

The methods described here are available from selected VARs. Value added reseller.



μSPE clean-up of QuEChERS extracts for GC and LC



• The μSPE cleanup achieves high quality results for diverse type of

analytes and foods (apple, kiwi, carrot, kale, orange, black olive,

The PAL syringe works as LC pump. Precisely

controlled flow rates in the load and elution steps


Lehotay et. al. Chromatographia 79, 17, pp 1113–30, 2016 http://link.springer.com/article/10.1007/s10337-016-3116-y/fulltext.html Morris, Schriner. J Agric Food Chem, 63, 5107-19, 2015 https://www.ncbi.nlm.nih.gov/pubmed/25702899


Fully integrated automated Bligh and Dyer extraction and dual-column analysis for metabolomics analyses of tissues and cells.


Emmanuel Varesio, Sandra Jahn, Sandrine Cudré, Gérard Hopfgartner, Life Sciences Mass Spectrometry, School of Pharmaceutical Sciences, University of Geneva, University of Lausanne, Switzerland; Renzo Picenoni, Guenter Boehm, Director Applications and Customer Communications, CTC Analytics AG, Zwingen, Switzerland

Conclusions

As laboratories are striving to uncover more «unknowns» and increase our understanding of biological processes there is a drive for procedures to become more efficient and repeatable. This is also true for extraction procedures which when performed manually can be time intensive and cumbersome, taking the valuable time of scientists. The automated Bligh and Dyer extraction described here was found to not only be more time efficient, but also to improve repeatability and quality of extraction and separation when compared to the standard manual approach.









Real time monitoring for off-flavour compounds using a functionalized autosampler with SPME-GC-MS/MS.


Jaewon Choi, Sung-Yun Ahn, Yuns Kim, Ilhwan Choi, Water Analysis & Research Center, K-water Wonkyoung Lee*, Moondon Choi*, Jongsu Park*, * Euro Science, Seoul, South Korea

Conclusions

- This real time monitoring system has been operating continously for several months. A temporary trend was observed for geosmin during 2 months.
- This system uses standard instrumentation, and is harmonized with the accredited method for drinking water of Korean Ministry of Environment (MOE) including adding internal standard and salt.







- 1 SPME Conditioning Station
- 2 10 μL Syringe Tool
- 3 10 mL Svringe Tool
- 4 Large Wash Module
- 5 Vial Tray Plate
- 6 Agitator Module
- 7 Vortex Mixer Module
- 8 Real Time Water Tank


Compound	Geosmin	2-MIB
R ²	0.998	0.995
Spiking conc. (ng/L)	1	2
LOD (ng/L)	0.16	0.17
LOQ (ng/L)	0.51	0.56
Precision (RSD %)	4.5	2.7
Recovery (%)	115	106

PAL Smart Consumables for Full Traceability

34

A fast, safe and reliable sample preparation is the key factor for high Full traceability productivity and reduces costs per sample: PAL3 Series II with Smart Smart Consumables status Lower cost of ownership Status display on a PAL Terminal Technology in combination with Smart Consumables provides the Higher sample throughput required process safety and efficiency. Highest result confidence **Smart Consumables** Chip on a tool **21**LSYSTEM Find more information about Smart Consumables





Error free usage (Plug and Play)

Full traceability

Tool

- Optimized Instrument up-time
- Status view
- Lower cost per sample
- Higher result confidence



Color code for easy identification of the syringe volume:




Find more information about **Smart Syringes**





PAL SMART SPME Arrow


The PAL SPME Arrow is the new patented technology for micro-extraction, combining trace level sensitivity with high mechanical robustness.





Color Code for easy optical identification of coating type and thickness


Optimized geometry:

- Up to 15 x more sorption phase
- Up to 6 x bigger surface
- Excellent extraction properties
- Designed for headspace and immersion sampling










PAL SMART SPME Fiber

PAL Smart SPME Fibers have been developed and optimized for the most successful SPME sampler, the PAL System Autosampler.



Color Code for easy optical identification of coating type and thickness



Integrated holder design facilitates installation:

- Automatic identification of SPME Fiber Type
- Excellent extraction properties
- Full traceability









PAL Accessories and Consumables - For safe and reliable processing.

The complete portfolio of PAL3 tested and certified consumables ensuring highest performance of every PAL System. PAL System Accessories are an integral part of the superior quality of every PAL System guaranteeing a safe and reliable operation.

- μSPE Cartridges for SPE of small volumes without the need for evaporation.
- PAL Vials and Caps for centrifugation, de- and recapping.
- PAL Pipette Tips for consistent pipetting with the PAL Pipette Tool.
- PAL Needle Seals for optimized LC injection with long lifetime, easy handling and no carryover.













Detailed information and specifications are listed in our Accessories and Consumables brochure.



Dimensions

	compact x-axis length	standard x-axis length	xtended x-axis length	xpanded x-axis length	xxpanded-axis length
	PALbase 534 X/Y unit	PALbase 850 X/Y unit	PALbase 1200 X/Y unit	PALbase 1600 X/Y unit	PALbase 2000 X/Y unit
Working Space	Width: 420 mm (16.8 inches) Depth: 255 mm (10.0 inches) Height: 420 mm (16.8 inches)	Width: 735 mm (28.9 inches) Depth: 255 mm (10.0 inches) Height: 420 mm (16.8 inches)	Width: 1090 mm (43.0 inches) Depth: 255 mm (10.0 inches) Height: 420 mm (16.8 inches)	Width: 1511 mm (60.4 inches) Depth: 255 mm (10.0 inches) Height: 420 mm (16.8 inches)	Width: 1915 mm (73.4inches) Depth: 255 mm (10.0 inches) Height: 420 mm (16.8 inches)
Footprint Instrument dimen- sions with standard legs	Width: 600 mm (24.0 inches) Depth: 795 mm (31.8 inches) Height: 770 mm (30.8 inches)	Width: 915 mm (36.6 inches) Depth: 795 mm (31.8 inches) Height: 770 mm (30.8 inches)	Width: 1270 mm (50.8 inches) Depth: 795 mm (31.8 inches) Height: 770 mm (30.8 inches)	Width: 1768 mm (70.7 inches) Depth: 795 mm (31.8 inches) Height: 770 mm (30.8 inches)	Width: 2172 mm (85.5inches) Depth: 795 mm (31.8 inches) Height: 770 mm (30.8 inches)
Sample Capacity	2 Tray Holders Up to 420 1 mL vials 324 2 mL vials 90 10/20 mL vials 6 MT/DW plates	4 Tray Holders Up to 840 1 mL vials 648 2 mL vials 180 10/20 mL vials 12 MT/DW plates	6 Tray Holders Up to 1260 1 mL vials 972 2 mL vials 270 10/20 mL vials 18 MT/DW plates	9 Tray Holders Up to 1890 1 mL vials 1458 2 mL vials 405 10/20 mL vials 27 MT/DW plates	13 Tray Holders Up to 2730 1 mL vials 2106 2 mL vials 58510/20 mL vials 39 MT/DW plates
Sample Capacity Temperature Con- trolled (4-40°C)	1 Peltier Stack Up to 420 1 mL vials 324 2 mL vials 90 10 mL vials 6 MT/DW plates	Max. 3 Peltier Stacks Up to 1260 1 mL vials 972 2 mL vials 270 10 mL vials 18 MT/DW plates	Max. 4 Peltier Stacks Up to 1680 1 mL vials 1296 2 mL vials 360 10 mL vials 24 MT/DW plates 9216 samples (with 54 x 384 well MTPs)	Max. 6 Peltier Stacks Up to 2520 1 mL vials 1944 2 mL vials 540 10 mL vials 36 MT/DW plates 13824 samples (with 54 x 384 well MTPs) Max. 27648 samples with 6 x Peltier Stack 12MT and 384 MTPs.	Max.9 Peltier Stacks Up to 37801 mL vials 29162 mL vials 810 10 mL vials 54 MT/DW plates 20796 samples (with 54 x 384 well MTPs) Max. 27648 samples with 6 x Peltier Stack 12MT and 384 MTPs.





Supported injection techniques

Tool	Ranges/Tools
Liquid Injection Tools	0.5 - 100 μ L syringes, 57 and 85 mm needle lengths 250 - 1000 μ L syringes, 57 and 85 mm needle lengths 5 & 10 mL syringes
uSPE Injection Tool	Liquid Tool for handling μSPE cartridges, holds syringes from 100-1000 μL
LCMS Tool	Injection volume LCMS-100, 1-80 μL Injection volume LCMS-250, 1-230 μL
Headspace Extraction	1000 μL / 2500 μL / 5000 μL, with corresponding tools
Solid Phase Micro Extraction (SPME & SPME Arrow)	SPME tool (holds commercially available fibers), SPME Arrow tool
Multiple Headspace Extraction (MHE)	MHE tool for 10 and 20 mL vials
ITEX Dynamic Headspace Extraction	ITEX DHS tool with Tenax TA adsorbent (other materials upon request)
Pipette Tool Injection	Direct injection from pipette tip into LC-valve

Available Modules

Agitator Module	Temperature controlled agitation, 40 - 200 °C, 250-750 rpm	
Barcode Reader Module	Reads 1D barcodes on 2, 10, 20 mL vials	
Centrifuge Module	Essential for sample prep, Combi, 2mL and 10 mL modules	
DeCapper Module	Opens & closes 2, 10 & 20 mL screw cap vials	
Dilutor Module	Fast and accurate dispensing of larger volumes of up to 5 different liquids	
Fast Wash Module	Efficient cleaning of syringes / needles (gauge 19-26) with 2 different solvents	
Fast Wash HF Module	Wash Module for flow rates up to 40 mL/min	
Flow Cell Module	Up to 6 flow cells for online sampling	
Heatex Stirrer Module	For mixing and heating in sample prep and SPME Arrow	
Large Wash Module	2 x 100 mL solvent container and waste port for cleaning of syringes / needles	
Liquid Cooler Module	Trayholder & tray for 32 10/20 mL vials (external cooler not included)	
Multiple Headspace Extraction Module	Tool & holder for MHE	
Peltier Stack Modules (2DW, 6DW, 12MT)	Temperature controlled storage 4 - 40 °C for a range vials and multititerplates	
Tray Cooler Module	Temperature controlled storage 4 - 40 °C for a range vials and multititerplates	
Solvent Module	3 x 100 mL solvent container for the addition of larger volumes of liquids	
SPME Arrow Conditioning Module	Conditioning of SPME Arrows and SPME Fibers, up to 350 °C, optional purge gas connection	
Standard Wash Module	4 x 10 mL solvent vial, 1 x 10 mL waste vial	
Valve Drive Module	For Rheodyne or Valco injection & switching valves	
Vortex Mixer Module	Efficient mixing for 2, 10, 20 mL vials	

LC Application Specifications

Туре	Specifications	Comment
Injection volume	0.1-10000 μL	Depending on syringe
Liquid injection, repeatability (gravimetrical)	Full loop < 0.1 % RSD Partial loop < 0.15 % RSD	20 μL loop, 4 x overfill 10 μL in 20 μL loop
Liquid injection, linearity (gravimetrical)	R > 0.9999	20, 40, 60, 80, 100 μL, n=3 each level
Liquid injection from small sample volume	3 injections 1 μL out of 5 μL	With bottom sense option and conical 150 μL vial
Carryover	< 30 ppm	With LC/MS Tool, blank measured after injection of CI-Hexidine 0.6 mg/mL

GC Application Specifications

Туре	Specifications	Comment
Injection volume	0.1-10000 μL	Depending on Syringe
GC liquid injection, repeatability	< 0.60 % RSD	Alkanes C14, C15, C16, 1 μL, split mode
GC liquid injection, linearity (gravimetrical)	R > 0.9999	20, 40, 60, 80, 100 μL, n=3 each level
GC head space injection, repeatability	< 1.00 % RSD	lso-octan, 10 μL in 20 mL vial, 500 μL injection
GC needle discrimination	C40/C20 > 0.98	Restek Florida Mix 1 µL, 100 ms fast split/splitless injection

47





Contact the experts for sample preparation:

BGB GC LC

Germany

www.bgb-shop.com

Switzerland BGB Analytik AG • Rohrmattstrasse 4 • 4461 Böckten • Phone +41 61 991 00 46 • Fax +41 61 991 00 25 • sales@bgb-analytik.com

Benelux BGB Analytik Benelux B.V. • Drielandendreef 42 • 3845 CA Harderwijk • Phone +31 341 700270 • Fax +31 341 700271 • sales.benelux@bgb-analytik.com

France BGB Analytik France S.A.S. • 81 Vie de l'Etraz • 01630 St. Jean de Gonville • Phone +33 450 488567 • Fax +33 450 562378 • sales.fr@bgb-analytik.com

BGB Analytik Vertrieb GmbH • Mühlestraße 1 • 79539 Lörrach • Phone +49 7621 5884270 • Fax +49 7621 5884289 • sales.de@bgb-analytik.com

For more information on the PAL RTC and RSI, including the latest application notes visit:

www.palsystem.com

