SunShell Core shell particle Sunniest Fully porous particle Bipheny


UHPLC & HPLC column

ChromaNik Technologies Inc.




Biphe<mark>nyl</mark>

Trifuctional silyl-reagent used



Feature of end-capping of Sunniest Biphenyl and SunShell Biphenyl (Sunniest End-capping)



Feature of Biphenyl

- The two benzene rings of the biphenyl group are twisted about 45 degree, and the biphenyl group shows different separation from the phenyl group alone.
- ✓ Even when an aqueous mobile phase containing no organic solvent is used, there is no decrease in retention, and separation with high reproducibility of retention time is achieved.
- \checkmark Greater retention of polar compounds such as caffeine.

Specification of Sunniest Biphenyl

ſ		Fully porous silica				Bonded phase				
		Particle size (µm)	Pore size (nm)	Specific surface area (m²/g)	Carbon loading (%)	Stationary phase	USP L line	End-capping	pH range	
	Sunniest Biphenyl	5	12	340	11	Biphenyl	L11	Sunniest end-capping	1.5 - 9	

Specification of SunShell Biphenyl

I		Core shell silica						Bonded phase				
		Particle size (µm)	Core size (µm)	Pore size (nm)	Specific surface area (m²/g)	Carbon loading (%)	Stationary phase	USP L line	End-capping	Maximum pressure	pH range	
	SunShell Biphenyl	2.6	1.6	9	150	5	Biphenyl	L11	Sunniest end-capping	60 MPa	1.5 - 9	

>Low bleeding

➢ High stability


Specification of other Biphenyl (cited from a brochure)

				-				-	
Company A Biphenyl	2.6	 10	100	11	Biphenyl	L11	TMS	60 MPa	1.5 - 8.5
Company B Biphenyl	2.7	 9	150	7	Biphenyldimethylsilane	L11	Yes	60 MPa	1.5 - 8.0
Company C Biphenyl	2.7	 9	150	7	Biphenyldimethylsilane	L11	Yes	60 MPa	1.5 - 9.0

1. Kinetex Biphenyl, 2.6 μm


Comparison of core shell Biphenyl columns Columns for comparison

Comparison of core shell Biphenyl phases using standard samples



Compared with other company's core shell type Bipjhenyl. Biphenyls from companies A, B and C showed similar separation patterns. It was confirmed that SunShell Biphenyl has a higher retention of amylbenzene No. 5 and a larger separation factor (hydrophobicity in the table) for the difference of one carbon between butylbenzene and amylbenzene, and is more hydrophobic than other company Biphenyl. The carbon loading is 5% for SunShell Biphenyl and 7% or more for all other Biphenyls, and the high hydrophobicity of SunShell, which has the lowest carbon content, indicates high density end-capping. Moreover, the separation of standard samples is very different for biphenyl and C18.

Comparison using pyridine



Pyridine is a compound that tends to tail due to residual silanol groups. When methanol is used as the organic solvent in the mobile phase rather than acetonitrile, tailing of pyridine is more likely to occur and the difference between brands becomes clear.



Biphenyl

Biphenyl

SunShell

Biphenvl

10

8


Company C

6

Retention time/min

Biphenvl

Comparison using formic acid



8-Quinolinol (oxine) is a metal chelating compound, and if there are metal impurities in the packing material, its peak becomes a tailing peak.

2

Flow rate: 1.0 mL/min

Detection: UV@250nm

Temperature: 40 ºC

4

Column dimension: 150 x 4.6 mm

Sample: 1 = 8-Quinolinol (Oxine)

2 = Caffeine


Mobile phase: CH₃CN/20mM H₃PO₄=10/90

n

Although most acidic compounds can be eluted without problems, formic acid is a typical acidic compound that easily causes tailing. Formic acid is tailed on Company A Biphenyl.




Comparison of core shell Biphenyl phases using basic compounds




Tailing of basic compounds (amitriptyline) occurs more easily in the acetonitrile/buffer mobile phase than in the methanol/buffer mobile phase. A mixture of acetonitrile and 20 mM phosphate buffer and a mixture of acetonitrile and 10 mM ammonium acetate were used as mobile phase to compare with core shell Biphenyls. Biphenyl from other companies had a terrible tailing for basic compounds, and it was also confirmed that the retention time was increased, which seems to be influenced by the adsorption to residual silanol groups.

Comparison of stability




Stability under acidic condition compared baseline changes with gradient elution of 1% aqueous phosphate solution and acetonitrile. Biphenyl groups are cut off from the silica surface under acidic condition and elute out of the column as the amount of organic solvent increases. Baseline variability detected cut Biphenyl groups, with SunShell showing the least desorbed and high acid stability. In addition, under weakly alkaline condition, the silica dissolved and the column-in side was dented, so the theoretical plate number of the columns were compared. Since the Biphenyl column has high reproducibility of retention time even in a 100% aqueous mobile phase (see page 6) and is effective for separating highly polar compounds, stability comparison was performed under the condition of pH 8 that does not contain an organic solvent. SunShell is more than twice as durable as other companies' columns.

Comparison of stationary phases using isomers of methylhippuric acid




Column: SunShell C18 2.6 μm 150 x 4.6 mm SunShell C30 2.6 μm 150 x 2.1 mm SunShell Phenyl 2.6 μm 150 x 4.6 mm SunShell Biphenyl 2.6 μm 150 x 4.6 mm Mobile phase: Organic solvent/25 mM Phosphate buffer pH 3.0 Flow rate: 1.0 mL/min, 0.2 mL/min for only C30



ChromaNik Techno


Temperature: 40 °C Detection: UV@230 nm Sample: *o-, m-, p*-Methylhippuric acid



The separation of ortho, meta and para isomers of methylhippuric acid was compared. The stationary phase was C18, C30, Phenyl-hexyl and Biphenyl, and the organic solvent in the mobile phase was methanol, acetonitrile, 2-propanol and tetrahydrofuran. The separation of metamethylhippuric acid and para-methyl hippuric acid changed depending on the organic solvent, and the separation was improved in the order of acetonitrile, methanol, 2-propanol, and tetrahydrofuran. Biphenyl showed a particularly high retention when alcohol was used as the organic solvent in the mobile phase, and the elution order of meta and para was reversed from that of other stationary phases. It is considered that this is due to the high hydrogen bonding capacity obtained on comparison of separation of standard samples. When acetonitrile is used, the π - π interaction between the solute and the stationary phase is weakened by the triple bond of CN in acetonitrile, so it is considered that the characteristics of Biphenyl cannot be fully exhibited. When tetrahydrofuran is used, tetrahydrofuran enters the stationary phase and a mixture of biphenyl group and tetrahydrofuran works as a stationary phase, so it is considered that the separation behavior was different from that when alcohol was used.


ChromaNik Technologies Inc

Comparison of stationary phases using creatinine, mandelic acid, hippuric acid and methylhippuric acid



Phosphate buffer was used as the mobile phase for separation of methylhippuric acid isomers on the previous page, but in this comparison of the separation including creatinine, the ammonium acetate buffer mobile phase that can be applied to LC / MS and the mobile phase in which an ion pair reagent were added to a phosphate buffer solution was used. 2-Propanol was used as the organic solvent. Under both mobile phase conditions, Biphenyl showed the highest retention and the best separation. A characteristic separation of Biphenyl is that the elution order of the isomers meta-methylhippuric acid and para-methylhippuric acid is reversed to that of C18, C8 and Phenyl-hexyl.

ChromaNik Technologies Inc.




Comparison of stationary phases using nucleic acid bases

Separation of nucleic acid bases was performed on C18, RP-AQUA and Biphenyl columns. Two mobile phases were used: methanol/10 mM ammonium acetate aqueous solution. Nucleic acid bases can be separated even in the mobile phase containing 20% methanol, but the retention time is shorter. In comparison of columns used this 20% methanol mobile phase, the retention factor (k5) of thymine at the 5th peak was highest in Biphenyl, and the resolution of the 5th and 6th peaks of thymine and adenine (R5,6) was also highest in Biphenyl. When an aqueous mobile phase containing no organic solvent was used, the retention factor of thymine eluting at last on C18 and RP-AQUA were 4.53 and 4.47, respectively. Biphenyl showed the largest retention factor for the compound eluted at last. The selectivity of Biphenyl is very different from that of RP-AQUA and C18, such as the order of sample elution. It has also been confirmed that the retention time stability when using an aqueous mobile phase is as high for Biphenyl as for RP-AQUA. The retention factor on C18 decreased to 15% after stopping the flow for 1 hour, confirming low reproducibility in retention.

★ Biphenyl stationary phase has a high hydrogen-bonding capacity and is effective for retaining polar compounds.

Moreover, separation with high reproducibility of retention time is achieved even with 100% aqueous mobile phase.





Ordering information of SunShell Biphenyl

Packings	Inner diameter (mm)	1.0	2.1	3.0	4.6	USP category		
	Length (mm)	Catalog number	Catalog number	Catalog number	Catalog number			
	30		C86931	C86331	C86431			
	50		C86941	C86341	C86441			
SunShell Biphenyl, 2.6 µm	75		C86951	C86351	C86451	L11		
	100		C86961	C86361	C86461			
	150		C86971	C86371	C86471			

Ordering information of Sunniest Biphenyl

Packings	Inner diameter (mm)	2.0	3.0	4.6	10	20	USP category	
	Length (mm)	Catalog number						
	50	E83241	E83341	E83441				
Supplect Diphonul Fum	100	E83261	E83361	E83461			111	
Sunniest Biphenyl, 5 μm	150	E83271	E83371	E83471			L11	
	250	E83281	E83381	E83481	E83781	E83881		

www.bgb-shop.com

Switz	erland:	BGB Analytik AG • Rohrmattstrasse 4 • 4461 Böckten • Phone +41 61 991 00 46 • Fax +41 61 991 00 25 • sales@bgb-analytik.com
		BGB Analytik SA • Route de Pré-Bois 20 • 1215 Genève 15 • Phone +41 22 788 49 43 • Fax +41 22 788 49 45 • sales.fr@bgb-analytik.com
Benel	lux:	BGB Analytik Benelux B.V. • Drielandendreef 42-44 • 3845 CA Harderwijk • Phone +31 341 700270 • Fax +31 341 700271 • sales.benelux@bgb-analytik.com
Franc	ce:	BGB Analytik France S.A.S. • 81 Vie de l'Etraz • 01630 St. Jean de Gonville • Phone +33 450 488567 • Fax +33 450 562378 • sales.fr@bgb-analytik.com
Germ	any:	BGB Analytik Vertrieb GmbH • Mühlestraße 1 • 79539 Lörrach • Phone +49 7621 5884270 • Fax +49 7621 5884289 • sales.de@bgb-analytik.com