

BIOANALYSIS CATALOG

GL Sciences Inc.

FROM SAMPLE PREPARATION TO ANALYTICAL SEPARATION

Proteomics

Spin Column Trypsin Digestion	MonoSpin Trypsin
Market Leading Phosphopeptide Enrichment	Titansphere 5 μm Bulk Media Titansphere 10 μm Bulk Media Titansphere Phos-TiO 10 μm Bulk Media Titansphere Phos-TiO MP Kit Titansphere Phos-TiO Kit Titansphere Phos-TiO for Large Volume Sample
Desalting Phosphopeptide-Enriched Samples	GL-Tip SDB and GC
Efficient Fractionation of Peptide Samples	GL-Tip SCX and SDB-SCX
Identification of Peptides/Proteins	Meter Scale MonoCap HighResolution 2000 Series Other MonoCap Standard Length Columns
Electrospray Emitter for ESI-LC/MS	MonoSpray

Sample Preparation for LC/MS

Solid Phase Extraction Spin Columns	MonoSpin and MonoSpin Large Type Series
Protein Removal & Phospholipid Removal Plates	Fast Remover Series
Desalting Phosphopeptide-Enriched Samples	GL-Tip SDB and GC
Efficient Fractionation of Peptide Samples	GL-Tip SCX and SDB-SCX

Exosome

High Purity Exosome Isolation Columns

EVSecond

Antibody

Rapid Purification of Antibodies

MonoSpin ProA and ProG

DNA

Purification and Extraction of DNA

MonoFas Series

Proteins

Spin Column Trypsin Digestion

MonoSpin Trypsin



Peptides

Market Leading Phosphopeptide Enrichment

Titansphere 5 μm Bulk Media Titansphere 10 μm Bulk Media Titansphere Phos-TiO 10 μm Bulk Media <u>NEW</u> Titansphere Phos-TiO MP Kit Titansphere Phos-TiO Kit Titansphere Phos-TiO for Large Volume Sample



Phosphopeptide

Desalting Phosphopeptide-Enriched Samples

GL-Tip SDB and GC

Efficient Fractionation of Peptide Samples


GL-Tip SCX and SDB-SCX



LC/MS/MS

Identification of Peptides/Proteins

Meter Scale MonoCap HighResolution 2000 Series Other MonoCap Standard Length Columns



Phosphorylated Protein Research


Protein phosphorylation is recognized as a fundamental process which regulates cell differentiation, growth, and migration. Analyzing protein phosphorylation is complicated by the low concentration of any given phosphoprotein and any one time, and the relatively low ionization efficiency of phosphopeptides in MS analysis. Therefore, enrichment of phosphopeptides and the relative reduction of non-phosphorylated peptides is critical to accurate analysis of protein digests by LC/MS.

GL Sciences' Titanium Dioxide (TiO2 or Titania) products have emerged as the most effect means of phosphopeptide enrichment of protein digests prior to LC/MS analysis, replacing IMAC as the primary means of phosphopeptide sample pretreatment. Enrichment by titanium dioxide and IMAC, remain, however, complimentary techniques and are often used in combination to obtain optimal phosphopeptide analysis.

What Makes GL Sciences' Titanium Dioxide Products Unique and Superior?


Titanium Dioxide exists in three crystaline forms, known as rutile, anatase, and brookite. Rutile and Anatase forms are the most common and most useful for phosphopeptide enrichment, and the ratio of rutile form to anatase form has significant implications for applicability to enrichment of phosphopeptides. GL Sciences' manufacturing technique for it's phosphopeptide enrichment products produces a highly spherical bead with the optimum ratio of crystal forms of TiO2. The primary reasons the GL Sciences' Pho-TiO products show superior performance is a direct result of the unique formulation of our titanium dioxide beads.

Principal of Phosphopeptide Enrichment



Phosphate groups are preferentially adsorbed to the surface of titanium dioxide under acidic conditions and are eluted under basic condition. Non-phosphorylated acid peptides non-specifically bound to the TiO2 can be reduced by adding acid modifiers to the loading and/or wash buffers.

Basics of Phosphopeptide Analyses by LC/MS/MS



Titansphere™ Bulk Materials

Bulk Sorbent Materials for Purification & Enrichment of Phosphopeptides

While GL Sciences' Phos-TiO spin columns based enrichment products are useful for most sample pretreatment applications, some investigators require bulk titanium dioxide media for specialized applications. Our market leading Titansphere Phos-TiO Kit is now available in bulk media which is Titansphere Phos-TiO 10 µm bulk media, and is optimized for purifying and enriching more phosphopeptide.

Applications

Efficient Purification from HeLa Cell Lysate

The data at right shows the superior performance of Titansphere TiO using the HeLa Cell Lysate consisting mainly of non-phosphorylated peptides. Titansphere TiO shows exceptional selectivity - almost 90% of the bound peptides were phosphopeptides, and excellent capacity for total phosphopeptide binding. A competitive TiO product is shown, binding mainly non-phosphorylated peptides and a much lower total number of discreet phosphopeptide species.

Sample: HeLa Cell Lysate, Sample volume: 50 µg, Titansphere TiO volume: 1 mg

Compare Titansphere TiO with IMAC

The graph at right shows how Titansphere TiO compares to an IMAC enrichment using Arabidopsis cell extract. Titansphere TiO provides substantially higher total capacity and a much higher number of discreet phosphopeptides isolated.

Sample: Arabidopsis Cell Extract, Sample volume: 100 $\mu\text{g},$ Titansphere TiO volume: 1 mg


Specifications


Description	Titansphere TiO
Particle Size	5 μm, 10 μm
Particle Shape	Spherical
Adsorption Spot	Titanium Dioxide Crystal
Pore Size	100 Å (10 μm)
pH Range	2 ~ 12
Gravity	1.74

Ordering Information

Titansphere[™] Bulk Materials (Previous Version)

Description	Particle Size	Qty	Cat. No.
Titansphere TiO 5 µm, 500 mg	5 µm	1/pk	5020-75000
Titansphere TiO 10 µm, 500 mg	10 µm	1/pk	5020-75010





Titansphere[™] Phos-TiO Bulk Material

Bulk Sorbent Materials for Purification & Enrichment of Phosphopeptides

While GL Sciences' Phos-TiO spin tips based enrichment products are useful for most sample pretreatment applications, some investigators require bulk titanium dioxide media for specialized applications. Our market leading Titansphere Phos-TiO Kit is now available in bulk media which is Titansphere Phos-TiO 10 µm bulk media, and is optimized for purifying and enriching more phosphopeptide.

References

- Phosphopeptide enrichment by aliphatic hydroxy acid-modified metal oxide chromatography for nano-LC-MS/MS in proteomics applications, Sugiyama N, Masuda T, Shinoda K, Nakamura A, Tomita M, Ishihama Y., Mol Cell Proteomics. 2007 Jun;6(6):1103-9.
- Highly selective enrichment of phosphorylated peptides using titanium dioxide, Nature Protocols 1, 2006, 1929-1935
- 3. Global, in vivo, and site-specific phosphorylation dynamics in signaling networks, Olsen JV, Blagoev B, Gnad F, Macek B, Kumar C, Mortensen P, Mann M., Cell. 2006 Nov 3;127(3):635-48.
- Successive and selective release of phospholylated Peptides captured by Hydroxy Acid-Modified Metal Oxide Chromatography, Yutaka Kyono, Naoyuki Sugiyama, Koshi Imami, Masaru Tomita, and Yasushi Ishihama, J Proteome Res., 2008, 7(10), 4585-93
- 5. Extended coverage of singly and multiply phosphorylated peptides from a single Titanium Dioxide Microcolumn, Anal Chem., 2015, 87(20), 10213-21

Ordering Information

Titansphere[™] Phos-TiO Bulk Material (Same Bulk Media packed into Titansphere Phos-TiO Kit)

Description	Particle Size	Qty	Cat. No.
Titansphere Phos-TiO Bulk 10um, 500 mg	10 µm	1/pk	5010-21315

NEW Titansphere[™] Phos-TiO MP Kit

Efficiently Enrich **BOTH** Singly and Multiply Phosphorylated Peptides

GL Sciences' is known as the best manufacturer for the enrichment/purification of phosphopeptide, which our Titansphere Phos-TiO kit and bulk resins are widely used throughout the world in major cancer research institutes and proteomics core facilities.

The new Titansphere Phos-TiO MP Kit employs a new protocol in the HAMMOC method, which enables highly efficient and selective recovery of not only singly, but also for multiply phosphorylated peptides. Specifically, the new kit fractions the singly and multiply phosphorylated peptides separately, which prevents ion suppression in LC-MS/MS detection and delivers higher recovery of multiply phosphorylated peptides.

Features

- High Recovery of not only Singly, but also for Multiply Phosphorylated Peptides.
- All Operation is done using an Easy-To-Use Centrifuge.


Sample Loading Capacity

Description	Content		
Sample	Tyr(PO ₃ H ₂) - Angiotensin II		
Spin Tip Sorbent Mass/Tip Volume	1 mg/ 200 μL 3 mg/200 μL		
Sample Loading Capacity	1.2 µg	3.5 µg	

Titansphere Phos-TiO MP Kit Contents

Cat. No.	5010-21282	5010-21283
Titansphere Phos-TiO MP Kit	24 pcs	24 pcs
Titansphere Sorbent Mass/Tip Volume	1 mg / 200 μL	3 mg / 200 µL
Spin Tip Quantity	24 pcs (6 x 4 packs)	24 pcs (6 x 4 packs)
Waste Fluid Tube Quantity	24 pcs	24 pcs
Recovery Tube (2.0 mL) Quantity	24 pcs	48 pcs
Recovery Tube (1.5 mL) Quantity	48 pcs	48 pcs
Solution B (Lactic acid) Quantity	2 mL	2 mL
Instruction Manual	1/pk	1/pk

NEW Titansphere[™] Phos-TiO MP Kit




Comparison of Recovery of **BOTH** Singly and Multiply Phosphorylated Peptides

The above experiment was done using HeLa cells 100 µg to prove the new Titansphere Phos-TiO MP Kit show highly efficient recovery of both singly and multiply phosphorylated peptides compared to SIMAC (Sequential Elution from IMAC) method.

The SIMAC method is generally separated into two major protocols, which the initial enrichment and separation of mono- and multi-phosphorylated peptides uses an Immobilized Metal ion Affinity Chromatography and a subsequent enrichment of the mono-phosphorylated peptides using titanium dioxide chromatography. Finally, the two separated phosphopeptide fractions are then analyzed by LC-MS/MS.

In the new Titansphere Phos-TiO MP Kit, singly phosphorylated peptides were recovered more when using an acidic elution buffer, while the multiply phosphorylated peptides were recovered more when using a basic elution buffer. As proven above, in total, the new Titansphere Phos-TiO MP Kit recovered more of both singly and multiply phosphorylated peptides than the SIMAC method.

Typical Enrichment Protocol using Titansphere Phos-TiO MP Kit



- Buffer A, Buffer B, Methylphosphonic acid, Disodium hydrogen phosphate, Ammonia aqueous solution and 5% Pyrrolidine aqueous solution are not included in the kit due to their designation as hazardous materials for the purposes of air transportation. Therefore, the user must prepare these solutions in their lab, according to the procedure detailed in the instruction manual.
- Solution B is included in the kit, which is Lactic acid.
- The preparation procedure of Buffer A and B are as follows.

	2% TFA solution	1mL
Buffer A	Acetonitrile	4mL
	Total	5mL (Use 3mL for making Buffer B)
	Solution B	1mL
Buffer B	Solution B Buffer A	1mL 3mL

Titansphere[™] Phos-TiO MP Kit

- Centrifuge Adapter, 24 pcs/pk (Cat. No. 5010-21514) must be purchased once to use the Titansphere Phos-TiO MP Kit.
- This centrifuge adapter is reusable.
- Lactic acid is Solution B, which is already included in the kit, however, can be purchased separately for future requirements.

Description	Sorbent Mass/Tip Volume	Quantity	Cat. No.
Titananhara Dhaa TiQ MD Kit	1 mg / 200 μL	24 pcs	5010-21282
Titansphere Phos-TiO MP Kit —	3 mg / 200 μL	24 pcs	5010-21283
Centrifuge Adapter	-	24 pcs	5010-21514
Lactic acid for Titansphere Phos-TiO (This is Solution B included in the kit)	15 mL	1/pk	5010-21295

Phos-TiO MP Kit with Desalting Columns

• These special packages includes optimized desalting columns/spin tips to be used with Phos-TiO MP Kit.

Package Contents	Package Cat. No.
Titansphere Phos-TiO MP Kit, 1mg /200 µL, 24 pcs/pk	
(Cat. No. 5010-21282)	
MonoSpin Large type C18, 30 pcs/pk	5010-21272
(Cat. No. 7510-11320)	5010-21272
GL-Tip SDB, 96 pcs/pk	
(Cat. No. 7820-11200)	

Package Contents	Package Cat. No.
Titansphere Phos-TiO MP Kit, 3mg /200 µL, 24 pcs/pk	
(Cat. No. 5010-21283)	
MonoSpin Large type C18, 30 pcs/pk	5010-21273
(Cat. No. 7510-11320)	5010-21275
GL-Tip SDB, 96 pcs/pk	
(Cat. No. 7820-11200)	

Titansphere Phos-TiO Spin Tips

• Spin tips are also available separately.

Description	Sorbent Mass/Tip Volume	Quantity	Cat. No.
Titansphere Phos-TiO Spin Tips	1 mg / 200	24 pcs	5010-21316
	1 mg / 200 μL	96 pcs	5010-21317
	2 mg / 200 l	24 pcs	5010-21307
	3 mg / 200 μL	96 pcs	5010-21308
9		•	

Titansphere[™] Phos-TiO Kit (Previous Version)

Enrichment of Phosphopeptide Using Spin Tips

This is the previous phosphopeptide enrichment kit that GL Sciences introduced to the market, which became the most popular kit worldwide. The titanium dioxide particles contained in the spin tips (available in 1 mg/10 μ L and 3 mg/200 μ L sizes) is specially treated to maximize selectivity for phosphorylated species, and the conditioning and washing buffers contain components to displace the few non-phosphorylated compounds which might originally adhere to the media.

Features


- High Recovery of Singly Phosphorylated Peptides.
- All Operation is done using an Easy-To-Use Centrifuge.

Sample Loading Capacity

Description	Content		
Sample	Tyr(PO ₃ H ₂) - Angiotensin II		
Spin Tip Sorbent Mass/Tip Volume	1 mg/ 10 µL	3 mg/200 µL	
Sample Loading Capacity	1.2 µg	3.5 µg	

Titansphere Phos-TiO Kit Contents

Cat. No.	5010-21309	5010-21310	5010-21311	5010-21312
Titansphere Phos-TiO Kit for Export	24 pcs	96 pcs	24 pcs	96 pcs
Titansphere Sorbent mass/Tip Volume	1 mg / 10 μL		3 mg / 200 μL	
Spin Tip Quantity	24 pcs (6 x 4 packs)	96 pcs (6 x 16 packs)	24 pcs (6 x 4 packs)	96 pcs (6 x 16 packs)
Waste Fluid Tube Quantity	24 pcs	96 pcs	24 pcs	96 pcs
Recovery Tube Quantity	Fube Quantity 24 pcs 96 pcs		24 pcs	96 pcs
Solution B (Lactic acid) Quantity	2 mL	6 mL	2 mL	6 mL
Instruction Manual	1/pk	1/pk	1/pk	1/pk




Phos-TiO Kits Outperform 4 Competitive TiO Based Products (MALDI-TOF/MS)

The data above show the purification efficiency of various TiO based products with a 2.5 µg sample of B-casein digest using MALDI-TOF/MS. Compared to the untreated condition (a), phosphopeptides were selectively purified when using Titansphere Phos-TiO Kit. Compared to competitive products (c to e) Titanpshere Phos-TiO Kit showed better selectivity. In general titanium dioxide is said that it has the worse adsorption efficiency of multi-phosphopeptides than IMAC. However, Titansphere Phos-TiO Kit showed higher selectivity, sensitivity and number of individual phosphopeptides isolated for 4 – phosphopeptides than IMAC (f). Metastable ion is a dephosphorylated peak.

Comparison between GL Sciences' Previous Version Bulk Media

Optimal TiO beads (Titansphere Phos-TiO Bulk 10um, Cat No. 5010-21315) are used for Titansphere Phos-TiO Kit. The existing Titansphere TiO bulk media were improved for better adsorption capacity of phosphopeptides. Compared to the existing Titansphere bulk media (ex: Cat No. 5020-75000, 5020-75010), Phos-TiO Kit showed 2.6 times more peak area and 1.6 times more identified phosphopeptides.



Ordering Information

Titansphere[™] Phos-TiO Kit

- Centrifuge Adapter, 24 pcs/pk (Cat. No. 5010-21514) must be purchased once to use the Titansphere Phos-TiO Kit.
- This centrifuge adapter is reusable.
- Lactic acid is Solution B, which is already included in the kit, however, can be purchased separately for future requirements.

Description	Sorbent Mass/Tip Volume	Quantity	Cat. No.
Titapaphara Dhac TiQ Kit	1 mg / 10 ul	24 pcs	5010-21309
Titansphere Phos-TiO Kit	1 mg / 10 μL	96 pcs	5010-21310
Titansphere Phos-TiO Kit	2 mg / 200 ml	24 pcs	5010-21311
	3 mg / 200 µL	96 pcs	5010-21312
Centrifuge Adapter	-	24 pcs	5010-21514
Lactic acid for Titansphere Phos-TiO (This is Solution B included in the kit)	15 mL	1/pk	5010-21295

Titansphere Phos-TiO Spin Tips

• Spin tips are also available separately.

Description	Sorbent Mass/Tip Volume	Quantity	Cat. No.
	1 mg / 10	24 pcs	5010-21302
Titopophoro Phoo TiO Coin Ting	1 mg / 10 μL	96 pcs	5010-21303
Titansphere Phos-TiO Spin Tips 3 mg / 200 μL	2 mg / 200	24 pcs	5010-21307
	96 pcs	5010-21308	

Titansphere[™] Phos-TiO for Large Volume Samples

Titansphere™ Phos-TiO for Large Volume Samples




Appropriate for Larger Scale/Volume Purifications

The same specialized bulk media used in our Phos-TiO Kit is available in 50 mg/3 mL and 100 mg/3 mL cartridges as an extension of the Phos-TiO product line.

Sample Loading Capacity

Description	Content		
Sample	Tyr(PO ₃ H ₂) - Angiotensin II		
Spin Tip Sorbent Mass/Tip Volume	50 mg/3 mL 100 mg/ 3 mL		
Sample Loading Capacity	50 µg	100 µg	

Typical Enrichment Protocol using Titansphere Phos-TiO for Large Volume Samples



Phos-TiO for Large Volume Samples Cartridges are intended for use with a desktop or other centrifuge. While some of the versions of Phos-TiO are resemble pipette tips or SPE cartridges, these products are not intended for use with pipettes or SPE vacuum manifolds; the cartridge internal configuration and particle size of the TiO bulk media requires centrifugal elution of all solutions.

Ordering Information

Titansphere[™] Phos-TiO for Large Volume Samples

Description	Sorbent Mass/Tip Volume	Qty	Cat. No.
Titansphere Phos-TiO	50 mg/3 mL	25 pcs	5010-21290
for Large Volume Samples	100 mg/3 mL	25 pcs	5010-21291

GL-Tip SDB and GL-Tip GC

Desalting Phosphopeptide-Enriched Samples Prior to LC-MS/MS

Phosphopeptides isolated using TiO2-based medias are typically desalted prior to analysis by LC-MS/MS, typically using a C18 (hydrophobic) micropipette tip. GL Sciences' SDB (styrene divinylbenzene) and GC (graphite carbon) centrifugeoperated micropipette GL-Tip retain more hydrophobic and hydrophilic peptides, respectively, than C18-based tips.


Features

GL-Tip SDB are more hydrophobic than C18 medias and allow retention of a wider range of phosphopeptides with high yield, allowing more accurate analysis of phosphopeptides species present in the sample. GL-Tip GC retain many more hydrophilic phosphopeptides than does C18; by using a combination of GL-Tip SDB and GC, almost all peptide samples can be desalted without sample losses due to lack of retention. Another highlight of this product is, the operation is very easy-to-use. Phosphopeptide-enriched samples are easily loaded, washed, and eluted using the same centrifuge-based technique used with Phos-TiO spin tips.

Sample Loading Capacity


Description	GL-Tip SDB	GL-Tip GC
Sample	Tyr(PO ₃ H ₂) - Angiotensin II	Gly-Gly-Tyr-Arg
Spin Tip Sorbent Mass/Tip Volume	200 µL	1 mg/200 µL
Sample Loading Capacity	60 µg	30 µg

Relative Retention of Peptides Collected using GL-Tip SDB and GC Desalting Tips



As illustrated above, the data indicating that GL-Tip SDB preferentially binds hydrophobic peptides while GC preferentially binds hydrophilic peptides.

Recommended Protocol using GL-Tip SDB and GL-Tip GC



GL-Tip SDB and GL-Tip GC

- Centrifuge Adapter, 24 pcs/pk (Cat. No. 5010-21514) must be purchased once to use the GL-Tip SDB and GL-Tip GC desalting spin tips.
- This centrifuge adapter is reusable.

Description	Tip Volume	Qty	Cat. No.
GL-Tip SDB	200 µL	96 pcs	7820-11200
GL-Tip GC	200 µL	96 pcs	7820-11201
Centrifuge Adapter	-	24 pcs	5010-21514





GL-Tip SCX and GL-Tip SDB-SCX


Spin Tips for Peptide Fractionation

GL-Tip SCX is packed with strong cation polymer (SCX) and GL-Tip SDB-SCX are packed with styrene divinylbenzene polymer (SDB) and strong cation polymer (SCX). GL-Tip SDB-SCX is packed in a two layer format consisting an SDB and SCX media. Undesalted peptide samples can be used in GL-Tip SDB-SCX as the first SDB layer can desalt the sample.

Sample Loading Capacity

Description	GL-Tip SCX	GL-Tip SDB-SCX
Sample	Angiotensin II	Angiotensin II
Tip Volume	200 µL	200 µL
Sample Loading Capacity	60 µg	60 µg


Comparison of Traditional Gradient Elution vs TFA Gradient Elution



A gradient elution using cation SCX media is commonly used in shotgun proteomics to fractionate peptide samples from complex samples such as cell or tissue extracts. The biggest challenge arises when identifying the same peptide from one fractionated peptide sample to another, which results in lowering efficiency. The newly developed TFA gradient elution method (patent applied) identifies more peptides without decreasing operation efficiency.

Comparison of Number of Quantified Peptides

Compared results between commercially available brand G's tip column. GL-Tip SCX recovered more peptides. The usage of the newly developed TFA gradient elution method provide less chance of identifying the same peptide from one fractionated sample peptide to another resulting in higher efficiency.



Number of Identified Non Redundant Peptides

Fractions	1	2	3	4	5	6	7	total
GL-Tip SCX	1996	1839	1684	1491	1311	1196	847	6085
Brand G	1552	1397	1004	1032	890	676	179	4704

(Data provided from the National Institute of Biomedical Innovation, Dr. Jun Adachi)

Ordering Information

GL-Tip SCX and GL-Tip SDB-SCX


- Centrifuge Adapter, 24 pcs/pk (Cat. No. 5010-21514) must be purchased once to use the GL-Tip SCX and GL-Tip SDB-SCX peptide fractionation spin tips.
- This centrifuge adapter is reusable.

Description	Tip Volume	Qty	Cat. No.
GL-Tip SCX	200 µL	96 pcs	7510-11203
GL-Tip SDB-SCX	200 µL	96 pcs	7510-11202
Centrifuge Adapter	-	24 pcs	5010-21514
Centrifuge Adapter for 96-Well Plate	-	1/pk	5010-21341
Centrifuge Adapter for 96-Well Plate	-	2 pcs	5010-21343



MonoCap HighResolution 2000 Series

Optimized for Identification of Peptides/Proteins for Proteome Research



MonoCap HighResolution 2000 is a 2 meter length monolithic silica capillary column which is designed for identifying extremely high number of peptides/proteins for proteome research via LC-MS/MS.

GL Sciences' MonoCap capillary columns, created synthetically via sol-gel method, and an octadecyl silane chemically bonded, has a very uniform three dimensional structure that shows excellent reproducibility from batch-to-batch. The solid structure of GL Sciences' monolithic silica eliminates the need for frits or filters at the ends of the column, thereby reducing dead volume that might otherwise lead to band broadening or sample recovery. The high porosity of our monolithic silica allows high flow rates to be used without loss of resolution or creation of high operating pressure. An optimized balance of through-pores and meso-pores provides the critically important combination of efficiency, separation speed, large volume sample-loading, and small volume sample-recovery.

MonoCap HighResolution provide extremely high efficiency, delivering over 200,000 plates for a 2,000 mm length column. The MonoCap HighResolution Ultra type deliver over 300,000 plates.

Physical Properties


Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 HighResolution 2000	Octadecyl Groups	15 nm	Yes	35 MPa
MonoCap C18 HighResolution Ultra 2000	Octadecyl Groups	11 nm	Yes	35 MPa
MonoCap HILIC-UP HighResolution 2000	Ureidopropyl Groups	12 nm	None	35 MPa

Based on monolithic technology, Merck KGaA, Darmstadt, Germany.



Discover New Peptides/Proteins

As proven below, MonoCap C18 HighResolution 2000 mm length column identifies simply more peptides/proteins compared to those traditional particle packed capillary HPLC columns.




_____20

For Identifying Highly Hydrophilic, Hydrophobic Peptides/Proteins

MonoCap HILIC-UP is an important addition to the MonoCap C18 HighResolution 2000 column series. MonoCap HILIC-UP can retain highly hydrophilic peptides/proteins which may lead to discovering new peptides/proteins where a C18 phase couldn't identify.

In HILIC, the higher the organic concentration, the greater the retention of more polar analytes. One of the biggest benefit of HILIC mode is, a high organic solvent concentration of the mobile phase will lead to a high sensitivity LC-MS/MS analysis.



for Proteomics LC-MS, K Horie et al. Anal. Chem. 2014, 86, 3817-3824

References

1. M.H.M. van de Meent et al.

Improvement of the liquid-chromatographic analysis of protein tryptic digests by the use of long-capillary monolithic columns with UV and MS detection, Anal Bioanal Chem, 2007,388, 195-200

2. Mio Iwasaki et al.

One-Dimensional Capillary Liquid Chromatographic Separation Coupled with Tandem Mass Spectorometry Unveils the Escherichia coli Proteome on a Microarray Scale, Anal. Chem. 2010, 82, 2616-2620

3. Mio Iwasaki et al.

Human Proteome analysis by using reversed phase monolithic silica capillary columns with enhanced sensitivity, J Chromatogr A 2012, 1228, 292-297

4. Ryota Yamana ea al.

Rapid and deep profiling of human induced pluripotent stem cell proteome by one-shot NanoLC-MS/MS analysis with meter-scale monolithic silica columns, J Proteome Res. 2013, 12, 214-21

5. Mari Ogawa-Ohnishi et al.

Identification of three hydroxyproline O-arabinosyltransferases in Arabidopsis thaliana, Nature Chem. Biol. 2013, 9, 726-730

- Satoru Okamoto et al. Root-derived CLE glycopeptides control nodulation by direct binding to HAR1 receptor kinase, Nature Commun. 2013,4, 2191
- 7. Kanta Horie et al.

Hydrophilic interaction chromatography using a meter-scale monolithic silica capillary column for proteomics LC-MS, Anal. Chem. 2014, 86, 3817-3824

MonoCap C18 HighResolution Ultra 2000

- End-fittings are not included.
- A column connection kit is available separately to ensure proper connections.
- Please refer to the below ordering information.

Description	I.D. (mm)	Length (mm)	Qty	Cat. No.
MonoCap C18 HighResolution Ultra 2000	10 µm	2000	1/pk	5020-10018

MonoCap C18 HighResolution 2000

- End-fittings are not included.
- A column connection kit is available separately to ensure proper connections.
- Please refer to the below ordering information.

Description	I.D. (mm)	Length (mm)	Qty	Cat. No.
MonoCap C18 HighResolution 2000	10 µm	2000	1/pk	5020-10015

MonoCap HILIC-UP HighResolution 2000

- End-fittings are not included.
- A column connection kit is available separately to ensure proper connections.
- Please refer to the below ordering information.

Description	I.D. (mm)	Length (mm)	Qty	Cat. No.
MonoCap HILIC-UP HighResolution 2000	10 µm	2000	1/pk	5020-10019

Connection Kit for MonoCap HighResolution 2000

- A dedicated connection kit for MonoCap C18 High Resolution 2000.
- Use this connection kit when connecting the column directly to the system.

Description		Cat. No.
1/16" PEEK Ferrule, SUS Nut, Sleeve, 2 pcs each. 1/32" PEEK Ferrule, SUS Nut, Sleeve, 2 pcs each.	1/pk	5020-10017



Zero Dead Volume Union

• Connect the tubing from the system to this union and install the column to achieve zero dead volume.

ю.	Cat. No.	Qty	Orifice Size	Description	
352	6010-72352	1/pk	0.25 mm		U-435
351	6010-72351	1/pk	178 µm		U-411
	001072	-/ pit	170 µm		0 111




23

MonoCap C18 Fast-Flow

Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure	
MonoCap C18 Fast-flow	Octadecyl Groups	15 nm	Yes	22 MPa	

Workable at a broad range of linear velocity from 0.5 to 5 mm/s without sacrificing efficiency and separation at high speed. The number of theoretical plates produced by MonoCap C18 Fast-Flow is nearly equivalent to a totally porous particle type capillary column packed with a 5 μ m packing material. Columns are protected by either metal or PEEK hardware.





MonoCap C18 Nano-flow



Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 Nano-flow	Octadecyl Groups	11 nm	Yes	22 MPa

MonoCap C18 Nano-flow produces higher number of theoretical plates compared to a totally porous particle type capillary column packed with a 3 μ m packing material. It can be operated at a wide range of flow rate with low back pressure and achieve very high sensitive results in Nano-LC-ESI/MS applications. Columns are protected by either metal or PEEK hardware.



MonoCap C18 WideBore Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure	
MonoCap C18 WideBore	Octadecyl Groups	11 nm	Yes	22 MPa	

The MonoCap C18 Fast-flow is also available in 0.5 mm I.D. size, which can be used at a wide range of flow rate from 6 to 100 μ L/min without sacrificing efficiency. The number of theoretical plates produced by MonoCap C18 WideBore is nearly equivalent to a totally porous particle type capillary column packed with a 5 μ m packing material. Columns are protected by a metal hardware.

MonoCap C18 Trap Column

Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 Trap Column	Octadecyl Groups	11 nm	Yes	20 MPa

MonoCap C18 Trap columns have a relatively big throughpore and workable at a high flow rate such as 10 μ L/min. This benefit makes MonoCap C18 Trap columns to be appropriate for on-line preconcentration or desalting of protein and peptide samples prior to HPLC separation with mass spectrometry detection. End-fittings are 1/16" (10-32 UNF). 1/32" end-fittings are also available upon request.



MonoCap Amide



Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure	
MonoCap Amide	Carbamoyl Groups	15 nm	None	22 MPa	


Amide groups are chemically bonded to the monolithic silica and makes it suitable for the analysis of sugars via HILIC mode. As the back pressure is significantly low, a 500 mm length MonoCap Amide column deliver over 40,000 plates offering high efficiency. Generally, HILIC mode uses acetonitrile at a concentration between 65-95 % in an aqueous buffer such as ammonium acetate or ammonium formate, which have high solubility in organic solvents. Columns are protected by either metal or PEEK hardware.

MonoCap SCX

Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap SCX	Benzenesulfonyl Groups	11 nm	None	20 MPa

MonoCap SCX is bonded with benzene sulfonic acid groups (strong cation exchange) and appropriate for 2D LC applications for the separation of biomolecules such as peptides and proteins.



MonoCap C18 Fast-Flow

- For end-fittings information, please refer to page 30.
- All 50 mm length PEEK columns does not come with a hardware and will be supplied with 3 pcs of columns only.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
	_	50		1/pk	5020-10102
		150	Metal	1/pk	5020-10101
	0.05	250		1/pk	5020-10100
	0.05 -	50	_	3/pk	5020-10002
		150	PEEK	1/pk	5020-10001
		250	_	1/pk	5020-10000
	_	50	_	1/pk	5020-10211
	_	150	Metal	1/pk	5020-10212
	0.075	250		1/pk	5020-10213
	0.075 -	50		3/pk	5020-10221
	_	150	PEEK	1/pk	5020-10222
MonoCap C18 Fast-Flow		250		1/pk	5020-10223
	_	50	_	1/pk	5020-10112
	_	150	Metal	1/pk	5020-10111
	0.1 -	250		1/pk	5020-10110
	0.1	50		3/pk	5020-10012
		150	PEEK	1/pk	5020-10011
	_	250		1/pk	5020-10010
	_	50	_	1/pk	5020-10122
		150	Metal	1/pk	5020-10121
	0.2	250		1/pk	5020-10120
	0.2 -	50		3/pk	5020-10022
	_	150	PEEK	1/pk	5020-10021
		250		1/pk	5020-10020

MonoCap C18 Nano-Flow

- For end-fittings information, please refer to page 30.
- All 50 mm length PEEK columns does not come with a hardware and will be supplied with 3 pcs of columns only.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
		50	– Motal –	1/pk	5020-10143
	0.05 -	150	- Metal -	1/pk	5020-10141
	0.05	50	– PEEK –	3/pk	5020-10043
		150	PEEN	1/pk	5020-10041
	_	50	- Motal -	1/pk	5020-10231
	0.075	150	- Metal -	1/pk	5020-10232
	0.075	50	– PEEK –	3/pk	5020-10241
ManaCan C18 Nana Flow		150	PEEK -	1/pk	5020-10242
MonoCap C18 Nano-Flow		50	Matal	1/pk	5020-10153
	0.1	150	- Metal -	1/pk	5020-10151
	0.1	50	PEEK -	3/pk	5020-10053
		150		1/pk	5020-10051
		50	- Metal -	1/pk	5020-10163
	0.2	150	Metal	1/pk	5020-10161
	0.2	50	– PEEK –	3/pk	5020-10063
		150	PEEN	1/pk	5020-10061

MonoCap C18 WideBore

• For end-fittings information, please refer to page 30.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
		50	Metal	_1/pk	5020-10202
MonoCap C18 WideBore	0.5	150	only	_1/pk	5020-10201
		250		1/pk	5020-10200

MonoCap C18 Trap Column

• For end-fittings information, please refer to page 30.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
		50	_	1/pk	5020-10026
		100	_ With Hardware	1/pk	5020-10038
	0.05	150		1/pk	NA
	0.05	50		1/pk	5020-10027
		100	_Without Hardware _	1/pk	5020-10039
		150		1/pk	NA
		50		1/pk	5020-10028
	0.075	100	_ With Hardware	1/pk	5020-10036
MonoCap C18 Trap Column		150		1/pk	NA
		50		1/pk	5020-10029
		100	_Without Hardware _	1/pk	5020-10037
		150		1/pk	NA
		50		1/pk	5020-10033
		100	_ With Hardware _	1/pk	NA
	0.2	150		1/pk	NA
	0.2	50		1/pk	5020-10034
		100	_Without Hardware _	1/pk	NA
			150		1/pk

MonoCap Amide

For end-fittings information, please refer to page 30. ٠

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
		150	Metal	1/pk	5020-10191
		250		1/pk	5020-10192
	0.075	500		1/pk	5020-10193
	0.075	150		1/pk	5020-10091
		250	PEEK	1/pk	5020-10092
		500		1/pk	5020-10093
	0.1	150		1/pk	5020-10181
		250	Metal	1/pk	5020-10182
MonoCap Amide		500		1/pk	5020-10183
Monocap Amide		150	PEEK	1/pk	5020-10081
		250		1/pk	5020-10082
		500		1/pk	5020-10083
		150		1/pk	5020-10171
	-	250	Metal	1/pk	5020-10172
	0.2	500		1/pk	5020-10173
	0.2	150	PEEK	1/pk	5020-10071
		250		1/pk	5020-10072
29		500		1/pk	5020-10073

MonoCap SCX

• For end-fittings information, please refer to the following information.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
		50		_1/pk	5020-10174
	_	150	Metal	1/pk	5020-10175
MonoCap SCX		250	Melai	1/pk	5020-10176
	0.2	500		1/pk	5020-10177
	0.2	50		1/pk	5020-10074
		150	חברע	1/pk	5020-10075
		250	PEEK	1/pk	5020-10076
		500		1/pk	5020-10077

End-fittings of MonoCap Monolithic Capillary HPLC Columns

Description	End-fittings Details
MonoCap C18 Fast-flow MonoCap C18 Nano-flow MonoCap C18 WideBore MonoCap Amide	 Metal Hardware Type End-fittings are Valco 1/16" (10-32 UNF). Valco 1/32" (6-40 UNF) end-fittings can also be arranged upon request, indicate "1/32" when ordering.
MonoCap SCX	 PEEK Hardware Type 1/16" male nut, ferrule and PTFE sleeve are included.

Connection Kit for MonoCap C18 Trap Column

Description	Cat. No.	
MonoCap C18 Trap Column Connection Kit 1/16" (Union Sleeve Capillary Tubing 2 pcs each, Nut Ferrule 4 pcs each)	5020-10044	
MonoCap C18 Trap Column Connection Kit 1/32" (Union Sleeve Capillary Tubing 2 pcs each, Nut Ferrule 4 pcs each)	5020-10045	
MonoCap C18 Trap Column Assembly Parts 1/16" (Nut Ferrule 4 pcs each)	5020-10046	
MonoCap C18 Trap Column Assembly Parts 1/32" (Nut [.] Ferrule 4 pcs each)	5020-10047	

MonoSpray

Monolithic Electrospray Emitter for ESI-LC/MS

MonoSpray is an electrospray emitter for ESI-LC/MS which a monolithic packing is packed into a fused silica sprayer offering numbers of benefits compared to those traditional sprayers packed with particle based packings. Frits are not installed in MonoSpray to keep the monolithic packing in place, which results in offering simply longer lifetime and avoiding bed splitting problems compared to those traditional sprayers packed with particle based packings. The very high porosity of monolithic packing allows a wide range of operational flow rates, even at high flow rates.

Features

- High Sensitivity Analysis
- High Chemical Stability
- High Physical Stability
- Wide Range of Operational Flow Rates

Ordering Information

MonoSpray FS

• For online Nano-ESI-LC/MS.

Description	Length (mm)	O.D. (µm)	I.D. (μm)	Qty	Optimum Flow Rate	Cat. No.
MonoSpray FS 50		370	F0 -	5 pcs	– ()]_]() /min –	5010-20001
			50	20 pcs		5010-20006
	50		75	5 pcs	- 0.2 - 2.0 μL/min -	5010-20002
			/5	20 pcs		5010-20007
			100	5 pcs		5010-20003
			100	20 pcs	– 1.0 - 5.0 μL/min -	5010-20008

Please inquire for other sizes.

MonoSpray C18 Nano

• Nano sprayer packed with octadecylated silica monolith offering reversed phased separation.

Description	Length (mm)	O.D. (µm)	I.D. (µm)	Qty	Optimum Flow Rate	Cat. No.	
MonoSpray C18 Nano	50	370	50 $\frac{1 \text{ pcs}}{4 \text{ pcs}}$	1 pcs	0.1 - 0.5 µL/min	5010-20011	
				4 pcs		5010-20016	
			75	1 pcs	0.1 - 1.0 µL/min	5010-20012	
				4 pcs		5010-20017	
				100	1 pcs		5010-20013
			100	4 pcs	0.3 - 2.0 μL/min	5010-20018	

Please inquire for other sizes.

Based on monolithic technology, Merck KGaA, Darmstadt, Germany.

MonoSpin Series


Low-Molecular Compounds Extraction and Purification

The low-pressure, high-flow, and low-liquid-retention properties of GL Sciences' monolith silica technology make it uniquely suited for handling of small samples. MonoSpin SPE centrifugal spin columns have been developed to improve concentration and yields in low-volume sample preparation.

Features


- Easy-to-Operate
- Ideal for Small Sample Volumes
- Wide Variety of Functional Groups
- Rapid Operation Time





How to Operate

Centrifuge elution allows loss-free and efficient processing of many samples simultaneously, with little or no liquid retained by the separation matrix. And, excellent mass transfer and rapid sample binding on MonoSpin's monolith silica allows extremely rapid sample preparation compared with other methods.



Formats

S Type (Small)

- Disc Size : φ 4.2 x 1.5 mm
- Sample Volume : 50 ~ 800 μ L
- Elution Volume : 50 ~ 800 µL
- Centrifugation Speed : 2,000 \sim 10,000 x g

L Type (Large)

- Disc Size : φ 9 x 3 mm
- Sample Volume : 0.5 ~ 8 mL
- Elution Volume : 0.5 ~ 8 mL
- Centrifugation Speed : $1,000 \times g$

Product Lineup

MonoSpin C18/C18 FF



Formats: S

Octadecyl functional group. Optimal for drug extraction in biological samples, and desalting & enrichment of peptide samples C18 FF type employs large through-pore monolith silica for high viscosity samples.


MonoSpin C18-AX



 $\begin{array}{c} \begin{array}{c} & & \\$

Bonded with octadecyl and trimethylaminopropyl, a mix mode type. Delivers great retention for high salt concentrated serum samples. Optimal for the recovery of acidic drugs.


MonoSpin SAX



Formats: S

Bonded with trimethylaminopropyl combining both strong anion exchange & weak hydrophobic interaction. Optimal for the extraction of acidic drugs.


MonoSpin NH2



Formats: S

Bonded with aminopropyl. Optimal for the enrichment of sugar chain and/or hydrophilic compounds by HILIC mode.

MonoSpin Amide




Formats: S

Bonded with amide. Optimal for the extraction of sugar chains and various hydrophilic acidic and basic compounds by HILIC mode.

MonoSpin TiO







Monolith skeleton coated with titanium dioxide. Excellent for the enrichment of phosphopeptides.

Product Lineup

MonoSpin ME





Bonded with iminodiacetic acid. Optimal for the recovery of trace metals.

MonoSpin Ph


Monolith

Formats: S

Phenyl functional group. Optimal for the recovery of hydrophobic drugs in biological samples due to its weak retentivity and different selectivity compared to a C18 phase.


MonoSpin C18-CX





Bonded with octadecyl and benzene sulfonic acid combining both ion exchange & hydrophobic interaction. Optimal for dissociated basic drug in biological samples. Delivers higher cleanup efficiency compared to C18 or SCX.


MonoSpin SCX



Formats: S

Bonded with benzenesulfonic acid combining both strong cation exchange & hydrophobic interaction. Optimal for the extraction of basic drugs.

MonoSpin CBA




Formats: S L

Bonded with carboxyl acid combining both weak cation exchange. Optimal for the extraction of basic drugs.

MonoSpin PBA






Specific column combined with phenyl boric acid. Excellent for the selective extraction of cis diol compounds, such as catechol amines.


Product Lineup

MonoSpin Trypsin




Immobilized trypsin is available for performing rapid and efficient tryptic digests of proteins.

MonoSpin Phospholipid



SL Monolith skeleton coated with TiO2 and ZrO2. Excellent for the adsorption and removal of phospholipids.

MonoSpin ProA



S 96 Formats:

Formats:

Formats:

S

Protein A immobilized affinity spin column for the rapid purification of antibodies.

MonoSpin ProG



Formats: S 96

Protein G immobilized affinity spin column for the rapid purification of antibodies.

S: Small Type L: Large Type

96: 96-well plate

Physical Properties

		S Type	(Small)	L Type (Large)		Curface		
Product	Stationary Phases	Through- pore (µm)	Meso- pore (nm)	Through- pore (µm)	Meso- pore (nm)	Surface Area (m ² /g)	Sample Loading Capacity (Small Type)	Filter
MonoSpin C18	Octadecyl	5	10	10	10	350	100 µg (Amitriptyline)	
MonoSpin C18 FF	Octadecyl	20	15	-	-	300	50 µg (Amitriptyline)	
MonoSpin Ph	Phenyl	5	10	-	-	350	100 µg (Amitriptyline)	
MonoSpin C18-AX	Octadecyl, Trimethylaminopropyl	5	10	-	-	350	100 µg (Ibuprofen)	
MonoSpin C18-CX	Octadecyl, Benzenesulfonic acid	5	10	-	-	350	100 µg (Amitriptyline)	
MonoSpin SAX	Trimethylaminopropyl	5	10	10	10	350	100 µg (Ibuprofen)	
MonoSpin SCX	Benzenesulfonic acid	5	10	10	10	350	100 µg (Amitriptyline)	
MonoSpin NH2	Aminopropyl	5	10	10	10	350	100 µg (Maltopentaose)	
MonoSpin CBA	Carboxyl	5	10	10	10	350	100 µg (Amitriptyline)	None
MonoSpin Amide	Amide	5	10	-	-	350	100 µg (Angiotensin Ⅱ)	
MonoSpin PBA	Phenyl boric acid	5	10	-	-	350	100 µg (Dopamine)	
MonoSpin TiO	Titanium dioxide	20	15	-	-	350	40 µg (Adenosine monophosphate)	
MonoSpin Trypsin	TPCK treated Trypsin	5	10	-	-	350	-	
MonoSpin ME	Iminodiacetic acid	5	10	10	10	350	25 µg (Cu ion)	
MonoSpin Phospholipid	TiO2 + ZrO2	5	10	10	10	350	10 µL (Human serum)	
MonoSpin ProA	Protein A	2	60	-	-	-	400 µg (Human IgG)	
MonoSpin ProG	Protein G	2	60	-	-	-	400 µg (Human IgG)	

Specifications

Description	MonoSpin S Type ^{*1}	MonoSpin FF ^{*2}	MonoSpin L Type
Disc Size	Φ 4.2 x 1.5 mm	Φ4.2 x 1.5 mm	Ф9 x 3 mm
Sample Volume	50 ~ 800 μL	50 ~ 800 μL	0.5 ~ 8 mL
Elution Volume	50 ~ 800 μL	50 ~ 800 μL	0.5 ~ 8 mL
Centrifugation Speed	2,000~10,000 x g	1,000 x <i>g</i>	1,000 x <i>g</i>
Sample Loading Capacity	100 µg	50 µg	1 mg



Appropriate for Various Viscosity Samples

MonoSpin series are ideal for the sample preparation of biological samples. MonoSpin C18 Fast Flow (FF) type is excellent for high viscosity biological samples. Select the appropriate MonoSpin column type depending on the viscosity of sample and volume.


L	ow	Sample Viscosity	High
MonoSpin C18 S Type (Small)	Ex: Urine	Ex: Serum, Plasma	Ex: Whole Blood
MonoSpin C18 FF Type (Fast Flow)			
MonoSpin C18 L Type (Large)			




MonoSpin Applications

The low-pressure, high-flow, and low-liquid-retention properties of GL Sciences' monolith silica technology make it uniquely suited for handling of small samples. MonoSpin SPE centrifugal spin columns have been developed to improve concentration and yields in low-volume sample preparation without requiring evaporation or reconstitution.


Purification and Enrichment of Trace Samples






MonoSpin Series

Purification of Amphetamine in Urine using MonoSpin C18




Recovery of Drugs in Serum using MonoSpin C18




Reproducibility study for drugs in serum sample for three days using MonoSpin C18 (n=10) MonoSpin demonstrated high reproducibility for purification of drugs.

	Sample	Con. (ng/mL)	Rec. (%)	RSD (%)		Sample	Con. (ng/mL)	Rec. (%)	RSD (%)		Sample	Con. (ng/mL)	Rec. (%)	RSD (%)
		5	91.2	4.8			5	83.7	3.9			5	83.7	7.0
	Desinversion	10	86.1	3.3		Paroxetine	10	84.1	7.8		Amitrintulino	10	81.8	2.8
	Desipramine	50	85.2	5.9			50	83.9	8.2		Amitriptyline	50	83.8	3.0
		250	88.4	6.5			250	86.7	7.5			250	88.4	2.7
		5	96.3	9.5			5	85.7	8.1			5	97.9	9.0
	Territoria	10	95.8	1.5		Maprotiline	10	84.7	3.2		Sulpiride	10	95.5	8.5
	Imipramine	50	94.5	0.9			50	88.6	5.4			50	90.8	2.6
		250	95.9	0.9			250	87.5	7.7			250	92.6	3.0
		5	96.8	11.6			5	106.3	9.9					
	El	10	87.1	5.0		Duloxetine	10	104.8	6.7					
	Fluvoxamine	50	86.8	8.1			50	99.8	8.7					
39 L		250	87.5	9.7			250	99.8	6.0					


Desalting of Protein Digests using MonoSpin C18



Recovery of Hormone in Serum using MonoSpin C18



Purification of Whole Blood using MonoSpin C18 FF (Fast Flow)



Features of MonoSpin C18 FF (Fast Flow)

MonoSpin C18 FF is ideal for high viscosity samples, such as whole blood and complex matrix samples.

Specification

Through-pore	20 µm	
Meso-pore	15 nm	
Disc Size	φ4.2 x 1.5 mm	
Sample Volume	50 ~800 μL	
Elution Volume	50 ~800 μL	
Centrifugation Speed	Under 1,000 x <i>g</i>	
Sample Loading Capacity	50 µg (Amitriptyline)	

MonoSpin C18 FF offer fast flow of viscosity samples at a low centrifugation speed $(1,000 \times g)$. The following is a comparison of flow of solvents between MonoSpin C18 and MonoSpin C18 FF.


Solvents	Volume	MonoSpin C18	MonoSpin C18 FF
Methanol	500 μL	0	0
Water	500 μL	400 µL	0
Serum*	500 μL	300 μL	0

Testing Conditions 1,000 x g 30 sec

* A supernatant from serum sample was used,

which was centrifuged at 10, 000 x g for 1 min.


Purification of Pyridylaminated (PA) Sugar Chain using MonoSpin NH2




Purified PA sugar chain by HILIC mode. MonoSpin NH2 additionally removes residual fluorescent labeling reagents.

Fractionation of Protein Digests using MonoSpin SCX


MonoSpin SCX provide a rapid and easy fractionation of peptides by stepwise elution using buffers with various salt concentration.




Purification of Paraquat and Diquat using MonoSpin CBA




Purification of Catecholamines using MonoSpin PBA



Purification of Organic Phosphorous Pesticides in Serum using MonoSpin TiO



Rapid Digestion of BSA using MonoSpin Trypsin



46

List of References


Products	Target Analytes	Sample Matrix	Concentration	Recovery Rate	Detection	Reference No.
	amitraz, metabolites	serum	5 ng/mL	95.5, 92.2 %	LC-MS	[1]
	dibudcaine, naphazoline	serum	5 - 10 ng/mL	70.2 - 78.6 %	LC-MS	[2]
	MA, AP, MDA, MDMA	urine	100 ng/mL	96 - 111 %	LC-UV	[3]
	9 cold medicines	serum	5 - 50 ng/mL	2.5 - 73.8 %	GC-MS	[4]
	amphetamines (AP, MA, MDA, MDMA)	urine	5 - 10 ng/mL	84 - 94 %	GC-MS	[5]
	eperison	serum	0.5 ng/mL	92.8 - 96.0 %	GC-MS	[6]
	paraquat, diquat, fenitrothion	serum, urine	25 - 100 ng/mL	51.3 - 106.1 %	GC-MS	[7]
Mana Cain® C10	arsenics	urine	1 ng/mL	91.9 - 106.5 %	GC-MS	[8]
MonoSpin [®] C18	MAM-2201	blood	1 ng/mL	-	LC-MS/MS	[9]
	a-PVP, a-PBP	urine	1 ng/mL	82 - 100 %	GC-MS	[10]
	a-PVP, a-PBP	hair	0.2 ng/mL	75.5 - 101.5 %	LC-MS	[11]
	Phthalic acid esters	physiological saline	0.2 - 50 μg/L	71.2 - 107.3 %	-	[12]
	<desalting></desalting>	digested peptides	-	-	-	[13]
	<desalting></desalting>	iTRAQ labeled samples	-	-	-	[14]
	MAM-2201	blood	2.5 - 100 ng/mL	1 ng/mL	-	[15]
	Naringin	grapefruit juice	10 - 500 µM	10 µM	-	[16]
	opiates benzodiazepines, metabolites	urine serum	10 ng/mL 1 - 10 ng/mL	69.2 - 98.9 % 83.3 - 112.3 %	LC-MS	[17]
MonoSpin [®] SCX	<pre-column derivatization="" fluorescence=""></pre-column>	-	-	-	-	[18]
	<desalting acid="" amino="" of=""></desalting>	-	-	-	-	[19]
ManaSpin® C18 CV	acidic and basic drugs	urine	1 - 25 ng/mL	65 - 123 %	GC-MS	[20]
MonoSpin [®] C18-CX	<halogenated compounds=""></halogenated>	cells	-	-	-	[21]
MonoSpin [®] C18-AX	amphetamines (AP, MA), opiates, THC	urine	2 - 10 ng/mL	93.1 - 108.1 %	GC-MS	[22]
MonoSpin [®] PBA	Adenosine	urine	6 µM	80 - 113 %	-	[23]

- [1] J. Chromatogr., B 867 (2008) 99-104.
- [2] J. Chromatogr., B 872 (2008) 186-190.
- [3] J. Chromatogr., A 1208 (2008) 71-75.
- [4] Chromatographia., 70 (2009) 519-526.
- [5] Anal. Chim. Acta., 661 (2010) 42-46.
- [6] J. Health Sci., 56 (2010) 598-605.
- [7] Anal. Bioanal. Chem., 400 (2011) 25-31.
- [8] J. Sep. Sci., 35 (2012) 2506-2513.
- [9] Forensic Toxicol., 31 (2013) 333-337.
- [10] Forensic Toxicol., 32 (2014) 68-74
- [11] J. Chromatogr., B 942-943 (2013) 15-20.
- [12] J Pharm Anal., 1 (2011) 92-99

- [13] Proteomics., 13 (2013) 751-755
- [14] Journal of proteomics., 84 (2013) 40-51
- [15] Forensic Toxicol., 31 (2013) 333-337
- [16] The Journal of Clinical Pharmacology., 54 (2013)
- [17] J. AOAC Int., 94 (2011) 765-774.
- [18] Biomed. Chromatogr., 26 (2012) 147-151.
- [19] Orig Life Evol Bjosph., 43 (2013) 99-108
- [20] J. Sep. Sci., 34 (2011) 2232-2239.
- [21] Toxicology., 314 (2013) 22-29
- [22] Forensic Toxicol., 31 (2013) 312-321.
- [23] Biosensors and Bioelectronics., 41 (2013) 379-385

Recovery of Metal Ions using MonoSpin ME


MonoSpin ME is bonded with iminodiacetic acid and optimal for the recovery and purification of metal ions. Specifically, it is excellent for the extraction and purification of trace Pb in blood or urine. Additionally, it is appropriate for removing inorganic divalent cations from sample to prevent ion suppression for LC-MS/MS applications.



Recovery rate of Cu²⁺ using Zeeman GF-A-AF system


Number of Injections	Volume of solvent introduced (mL)	Recovery rate (%)
1	0.8	98±4
2	1.6	97±5
3	2.4	95±5
4	3.2	95±5
5	4	94±3

Retention Characteristics of Metal Element using Iminodiacetic Acid Functional Groups with Various pH



Removal of Phospholipids using MonoSpin Phospholipid


MonoSpin Phospholipid removes more than 90% of phospholipids from biological samples resulting in eliminating ion suppression in LC-MS/MS analysis. The MonoSpin Phospholipid also removes phospholipids from a serum sample volume of 50 μ L.



Phospholipid Removal Efficiency of MonoSpin Phospholipid

Retention Mechanism of Phospholipids

Monolith skeletal structure coated with TiO_2 and ZrO_2 selectively interacts with metal oxides and phosphorylated compounds, resulting in removing more than 90 % of phospholipids.




Rapid Purification of Antibodies using MonoSpin ProA and ProG

MonoSpin ProA and MonoSpin ProG are immobilized with protein A or protein G onto a silica monolith offering rapid purification of antibodies. A 96-well plate format is also available for high throughput purification.

Features

The silica is modified with a hydrophilic polymer and then immobilized with either Protein A or Protein G to prevent the adsorption of proteins, resulting in rapid purification and high recovery of antibodies.




Specification				
Bonded Phase	Protein A or Protein G			
Through-pore	2 µm			
Meso-pore Size	60 nm			
Disc Size	Φ4.6 x 1.5 mm			
Sample Volume	50 - 500 μL			
Decovery Date	MonoSpin ProA: IgG 90 % (With 400 mg IgG)			
Recovery Rate	MonoSpin ProG: IgG 90 % (With 300 mg IgG)			
Elution Volume	50 μL			
Centrifugation speed	2,300 × g			

Antibody Compatibility Table


Species	Antibody Class	Protein A	Protein G
	IgG	Ø	Ø
	IgG1	Ø	Ø
	IgG2	Ø	Ø
	IgG3	_	Ø
	IgG4	Ø	Ø
Human	IgM	—	—
	IgA	—	—
	IgE	—	—
	IgD	—	—
	Fab	0	0
	ScFv	0	_






Comparison of Elution Volume & Recovery Rate with other Brands' Products

MonoSpin ProA only requires 100 μ L elution buffer, providing a recovery rate of at least 90% IgG. On the other hand, other brands' products requires 400 μ L of elution buffer with a recovery rate of 70% IgG.



Recovery Rate and Reproducibility of IgG from medium cultured CHO cells with MonoSpin ProA 96-Well Plate



Ordering Information

MonoSpin S Type (Small) Columns

• Each MonoSpin S Type (Small) columns are attached with 1.7 mL recovery tubes and 2.0 mL waste tubes.

Description	Qty	Cat. No.
ManaSpin C19	50 pcs	5010-21700
MonoSpin C18 —	100 pcs	5010-21701
ManaSpin C19 EE	50 pcs	5010-21670
MonoSpin C18 FF —	100 pcs	5010-21671
MonoSpin Dh —	50 pcs	5010-21733
MonoSpin Ph —	100 pcs	5010-21734
MonoSpin C19 AV	50 pcs	5010-21735
MonoSpin C18-AX —	100 pcs	5010-21736
MonoSpin C18-CX —	50 pcs	5010-21731
	100 pcs	5010-21732
MonoSpin SAV	50 pcs	5010-21720
MonoSpin SAX —	100 pcs	5010-21721
MonoSpin SCX —	50 pcs	5010-21725
	100 pcs	5010-21726
MonoSpin NH2 —	50 pcs	5010-21710
	100 pcs	5010-21711
MonoSpin CBA —	50 pcs	5010-21729
	100 pcs	5010-21730
MonoSpin Amide —	50 pcs	5010-21727
	100 pcs	5010-21728
MonoSpin PBA —	50 pcs	5010-21715
	100 pcs	5010-21716
MonoSpin TiO —	50 pcs	5010-21705
	100 pcs	5010-21706
* MonoSpin Trypsin —	50 pcs	7820-11300
	100 pcs	7820-11301
MonoSpin ME —	50 pcs	5010-21737
	100 pcs	5010-21738
MonoSpin Phospholipid —	50 pcs	5010-21698
	100 pcs	5010-21699

* MonoSpin Trypsin must be refrigerated when not in use.







MonoSpin S Type (Small)

Recovery Tube (1.7 mL) Waste Tube (2 mL)

Ordering Information

MonoSpin L Type (Large) Columns

Description	Qty	Cat. No.
MonoSpin L C18	30 pcs	7510-11320
MonoSpin L SAX	30 pcs	7510-11321
MonoSpin L SCX	30 pcs	7510-11322
MonoSpin L NH2	30 pcs	7510-11323
MonoSpin L CBA	30 pcs	7510-11324
MonoSpin L ME	30 pcs	7510-11325
MonoSpin L Phospholipid	30 pcs	7510-11326



* Each MonoSpin L Type (Large) columns does not come with any recovery and waste tubes.

* Please prepare a 50 mL centrifuge tube separately (Ex: Falcon tube).

MonoSpin ProA, MonoSpin ProG

Description	Qty	Cat. No.
MonoSpin ProA	10 pcs	7510-11310
MonoSpin ProG	10 pcs	7510-11311
MonoSpin ProA 96-Well Plate	1/pk	7510-11312
MonoSpin ProG 96-Well Plate	1/pk	7510-11313

* MonoSpin ProA, ProG must be refrigerated when not in use.

MonoSpin S Type (Small) Trial Kits

• The following trial kits are available for purchase to test a whole range of MonoSpin columns to make the best decision on which MonoSpin to use.

Description	Available Phases	Qty	Cat. No.
MonoSpin Trial Kit 1	C18, TiO, SCX, SAX, 10 pcs each.	10 pcs/4 pk	5010-21740
MonoSpin Trial Kit 2	C18, Amide ,CBA, NH2, 10 pcs each.	10 pcs/4 pk	5010-21741
MonoSpin Trial Kit 3	SCX, SAX, CBA, NH2, 10 pcs each.	10 pcs/4 pk	5010-21742

- MonoSpin Trial Kit 1: Optimal for drug extraction in biological samples & purification of pesticides.
- Monospin Trial Kit 2: Compatible with both hydrophilic/hydrophobic applications. Optimal for purification of peptide and sugar chains.
- MonoSpin Trial Kit 3: Optimal for purification of ionic analytes.

FastRemover for Protein

Maximizes Sample Yield


FastRemover is a 96-well type filter plate ideal for preparing precipitated protein samples. High-throughput processing of plasma samples is performed simply, accurately, and reproducibly.

Features

- Easy filtration of biological samples.
- Trace analytes can be processed with minimal sample loss owing to the low volume design of the elution tip and filter.
- Perfect for processing with automated vacuum instruments.
- High sensitivity analysis is unaffected by contamination from plasticizers or other impurities found in other 96-well plates.
- Removal of microparticle contaminants enables injection to LC/MS/MS directly from the collection plate.


Typical Protocol using FastRemover for Protein

To demonstrate the performance of FastRemover for Protein, a BSA solution was prepared as follows:



Adsorption Test

A standard mixture containing 7 compounds were analyzed to evaluate potential non-specific adsorption to the plate. As shown in the following chromatograms, FastRemover for Protein provides minimal loss of target samples.



Ordering Information

FastRemover for Protein

Description		Qty	Cat. No.
EastDomover for Drotein (0.45 um)	96-well	1/pk	7820-11001
FastRemover for Protein (0.45 µm)		5/pk	7820-11005
FastRemover for Protein (0.20 µm)	96-well	1/pk	7820-11011
		5/pk	7820-11015

Related Accessories

Description	Qty	Cat. No.
Vacuum Manifold with shims	1 Set	5010-33101
Sealing Mat for Microplate, WSM-3SX (PTFE/SILCON)	5/pk	1030-43831
Sealing Tape for Microplate, (Polyolefin)	100/pk	1065-70002

FastRemover for Phospholipid

Rapid and Efficient Removal of Proteins and Phospholipids


The FastRemover for Phospholipid 96-well plate deliver a rapid and effective removal of proteins and phospholipids in plasma and serum samples without sacrificing the recovery of your target analytes.

Features

- Simple & easy protocol to remove proteins and phospholipids.
- High sensitivity analysis is unaffected by contamination from plasticizers or other impurities found in other 96-well plates.
- Removal of microparticle contaminants enables injection to LC/MS/MS directly from the collection plate.
- Removes more than 90% of phospholipids resulting in eliminating ion-suppression.
- Prolong HPLC/UHPLC column lifetime by removing proteins and phospholipids that can damage your column.


Typical Protocol using FastRemover for Phospholipid

The presence of phospholipids in plasma or serum samples is one of the major problems in LC/MS-(MS) analysis. Phospholipids can build up on your MS system and bleed off the HPLC/UHPLC column, causing ion suppression, shifts in retention time and peak shape and necessitating time consuming column and system maintenance. Use of FastRemover for Phospholipid 96-well plate will eliminate these effects and extend the lifetime of your HPLC/UHPLC column and deliver more predictable/accurate mass spectrometry results.




Extend HPLC/UHPLC Column Lifetime

Over the course of multiple injections, phospholipids build up and can lead to reduced column lifetime, showing increase in column back pressure, decrease in column sensitivity and efficiency. The figure on the right illustrates the removal efficiency of phospholipids, proteins and microparticles by FastRemover for Phospholipid.




Industry Leading High Recovery for Bioanalysis

As shown below, the FastRemover for Phospholipid 96-well plate deliver a rapid and effective removal of proteins and phospholipids in plasma and serum samples without sacrificing the recovery of your target analytes.



Comparison of Recovery Rate using Various Solvents for Deproteinization


The following seven solvents were used to deproteinate a serum sample. As a result, 0.1% formic acid in 100% acetonitrile showed the best recovery for not only basic, but also for acidic compounds.



Ordering Information

FastRemover for Phospholipid

Description	Qty	Cat. No.
FastRemover for Phospholipid (0.2 µm)	1/pk	7510-11021



EVSecond

Exosome Purification Columns

Recent studies have reported significant roles of extracellular vesicle "Exosome" in development and progression of various diseases including cancer metastasis. Therefore exosomes are considered as attractive targets for biomarkers and drug development. However, it remains difficult to isolate high-purity exosomes from biological fluids such as serum. EVSecond is a size exclusion chromatography open column optimized for effective purification of exosomes. Highly-purified exosomes can be easily collected from serum, plasma, or cell culture supernatant.


Features

- Simple gravity-flow handling without ultracentrifugation.
- EVSecond-purified exosomes possess efficient purity for comprehensive miRNA, proteome, and metabolome analysis.
- Exosomes are gently eluted in PBS without structural damage, allowing re-administration experiments of collected exosomes to cells or animals.

Advantages Over Traditional Procedures

- Much higher-purity exosomes can be obtained compared to ultracentrifugation or polymer precipitation methods.
- Unlike immuno-affinity purification using anti-tetraspanin antibodies, whole exosomes can be collected regardless
 of surface antigen profiles.


Typical Protocol using EVSecond



61

Purification of Exosomes from Human Serum

A large amount of free proteins, metabolites, and nucleotides are involved in serum samples. Insufficient purification of exosomes often causes co-detection of non-exosomal components, leading to incorrect quantification results in omics studies. Exosomes were isolated from 100, 200, 300, or 700 µL of human serum using EVSecond method. Exosomes were clearly separated from serum free proteins such as albumin or immunoglobulins.



Ordering Information

EVSecond

Description	Qty	Cat. No.
	10 pcs	5010-21390
EVSecond —	25 pcs	5010-21392
GL-SPE EXO Fraction Rack	1 set	5010-50450

EVSecond was developed based on the cooperation from Dr. Koji Ueda from Graduate School of Frontier Sciences, The University of Tokyo.


MonoFas DNA Purification Kit I

DNA Extraction & Purification


MonoFas DNA Purification Kit I purifies DNA from PCR products and agarose gels. Purified DNA can be used for sequencing, cloning/ligation, restriction digests, etc.

Features


- Multiple Roles Purifies DNA from PCR matrices or agarose gels.
- Rapid purification in about 4 minutes.
- High recovery rates even from sample volumes as low as 10 µL.
- Purify DNA fragments from 35 bp to 35 kbp.



Purify DNA fragments from 35bp up to 35kbp



Large Quantities of Agarose Gels can be Processed




MonoFas routinely extracts DNA from up to 1 g of agarose gel at once.

400 bp M : pHY Marker 1 : Unpurified PCR products 2 : Purified from agarose gel (MonoFas DNA Purification Kit I)

Specifications

Description	Purification from PCR products	Extraction from agarose gel
Time	4 mins	9 mins
Maximum DNA Binding Amount	<10 µg	<10 µg
Maximum Agarose Gel Throughput	-	<1g
Minimum Elution Amount	10 µL	10 µL
Column Volume	1 mL	1 mL
Processable DNA Range	35 bp - 35 kbp	35 bp - 35 kbp
	>85 %(100 bp - 5 kbp)	>80 %(100 bp - 5 kbp)
Recovery Rate	>60 %(5 kbp - 35 kbp)	>50 %(5 kbp - 35 kbp)
Primer Removal Percentage	95 %	_

Typical Protocol using MonoFas DNA Purification Kit I





Condition: Cycle sequencing method with Big Dye Terminator v3.1 manufactured by ABI Model: ABI 3730 Genetic Analyzer



Constant Centrifuge Acceleration $> 2,000 \times g (19,600 \text{ m/s}^2)$

2 mins : DNA purification from PCR products 7 mins : DNA purification from agarose gel

Ordering Information

MonoFas DNA Purification Kit I

Description	Qty	Cat. No.
	50 pcs	5010-21530
MonoFas DNA Purification Kit1, EXPORT	100 pcs	5010-21531
	250 pcs	5010-21532




MonoFas DNA Purification Kit III

Plasmid DNA Extraction from E.coli

MonoFas Plasmid Extraction Kit III is designed to purify plasmid DNA from E. Col cultures. The extracted plasmid DNA can be used without further purification for sequence analysis, restriction digestion, cloning/ligation, etc.

Features

- Rapid Plasma Purification in 8 minutes.
- Highly Purified Plasmid DNA.
- Wide Variety of Functional Groups
- BAC Clone Purification can also be performed.
- Stable Recovery Rate



Ordering Information

MonoFas DNA Purification Kit III

Description	Qty	Cat. No.
	50 pcs	5010-21533
MonoFas DNA Purification Kit3, EXPORT	100 pcs	5010-21534
	250 pcs	5010-21535

MonoFas BAC Extraction Kit V

Efficient Isolation of BAC, Cosmid DNA

MonoFas Cultured Cell DNA Extraction Kit VI

Isolation of High Yields of Pure DNA in 5 mins without Desalting & Damaging DNA

MonoFas DNA Bacteria Extraction Kit VII

DNA Extraction from Gram Positive & Negative Bacteria

MonoFas DNA Buccal Swabs Extraction Kit VII

DNA Extraction from Buccal Swabs without Desalting

MonoFas DNA Mouse and Rat Tail Extraction Kit IX

Rapid Extraction of DNA from Mouse & Rat Tails

MonoFas DNA Stool Extraction Kit X

Stool Extraction Kit Designed to Extract Bacteria & Epithelial Tissues in Intestine

MonoFas DNA Processed Food Extraction Kit XI

Rapid Extraction of DNA from Processed Food

MonoFas Plant DNA Extraction Kit XII

Rapid Extraction of DNA from Processed Food

For products details, please inquire.

Based on monolithic technology, Merck KGaA, Darmstadt, Germany.

InertSearch

InertSearch Application Notes

Access to the latest pharmaceutical, life science, environmental and food applications at

www.glsciences.com/tech/inertsearch



Worldwide Ordering Information

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255 Torrance, CA 90503 Phone: 310-265-4424 Fax: 310-265-4425 Email: <u>info@glsciencesinc.com</u> Web: www.glsciencesinc.com

GL Sciences B.V.

De Sleutel 9 5652 AS Eindhoven The Netherlands Phone: +31 (0)40 254 95 31 Email: info@glsciences.eu Web: www.glsciences.eu

International Distributors

Visit our Website at www.glsciences.com/distributors

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome Shinjuku-ku, Tokyo, 163-1130, Japan Phone: +81-3-5323-6620 Fax: +81-3-5323-6621 Email: world@gls.co.jp Web: www.glsciences.com

The specification are subject to change without notice due to continual improvements.

www.bgb-shop.com

 Switzerland
 BGB Analytik AG • Rohrmattstrasse 4 • 4461 Böckten • Phone +41 61 991 00 46 • Fax +41 61 991 00 25 • sales@bgb-analytik.com

 Benelux
 BGB Analytik Benelux B.V. • Drielandendreef 42 • 3845 CA Harderwijk • Phone +31 341 700270 • Fax +31 341 700271 • sales.benelux@bgb-analytik.com

 France
 BGB Analytik France S.A.S. • 81 Vie de l'Etraz • 01630 St. Jean de Gonville • Phone +33 450 488567 • Fax +33 450 562378 • sales.fr@bgb-analytik.com

 Germany
 BGB Analytik Vertrieb GmbH • Mühlestraße 1 • 79539 Lörrach • Phone +49 7621 5884270 • Fax +49 7621 5884289 • sales.de@bgb-analytik.com