

HAMILTON

Hamilton HPLC Column Compound Index

**Save Time...
Look inside before
you buy an HPLC
column!**

**Find the compound
you need to separate
in the Compound
Index and request
a copy of the
application
chromatogram.**


HAMILTON

HPLC APPLICATIONS

Hamilton Company
P.O. Box 10030
Reno, Nevada 89520-0012
Toll-Free 800-648-5950
Tel: +1-775-858-3000
Fax: +1-775-856-7259
e-mail: sales@hamiltoncompany.com
<http://www.hamiltoncompany.com>

Application #132 Fluorescein
Labeled 12mer
150 x 4.1 mm PRP-1 (P/N 79425)
1 Failure Sequences
2 Fluorescein Labeled 12mer

Conditions: A) 20 mM Sodium Acetate pH 6.7;
B) Acetonitrile. Linear Gradient 0-50%
B in 20 min. Ambient.
Flow: 1.0 mL/min.
Injection: 5 μ L
Detection: UV at 260 nm



HAMILTON


HPLC APPLICATIONS

Hamilton Company
P.O. Box 10030
Reno, Nevada 89520-0012
Toll-Free 800-648-5950
Tel: +1-775-858-3000
Fax: +1-775-856-7259
e-mail: sales@hamiltoncompany.com
<http://www.hamiltoncompany.com>

Application 308
Demeclocycline by USP Method
250 X 4.1 mm PRP-1 (P/N 79422)

1. Demeclocycline, 1000 ppm

Conditions: 80g t-Butyl alcohol,
100 mL 0.2 M Phosphate buffer
(pH 9.0), 150 mL 0.02M Tetrabutyl-
ammonium hydrogen sulfate, 100 mL
0.01M EDTA diluted to 1L with DI water.
Isocratic. 60°C
Flow: 1 mL/min
Injection: 10 μ L
Detection: UV at 254 nm



To receive a copy of the application chromatograms of interest, e-mail: hplc@hamiltoncompany.com, or call us **800-648-5950 ext. 290** or complete the enclosed reply card.

A

1.	AAF-Modified Oligonucleotide 9mer
2.	AAF-Modified Oligonucleotide 11mer
3.	AAF-Modified Oligonucleotide 12mer
4.	Acepromazine
5.	Acetamide
6.	Acetaminophen
7.	Acetate
8.	Acetic Acid
9.	Acetoin
10.	Acetone
11.	Acetophenazole
12.	Acetophenetidine
13.	Acetophenone
14.	N-Acetyl-2-aminofluorene Modified 9mer
15.	N-Acetyl-2-aminofluorene Modified 11mer
16.	N-Acetyl-2-aminofluorene Modified 12mer
17.	Acetylcholine
18.	Acetylferrocene
19.	N-Acetyl-D-Galactosamine
20.	Acetyl-Glucosamine
21.	N-Acetyl-D-Glucosamine
22.	Acetylsalicylic Acid
23.	Acid Blue Dye 1
24.	Acid Blue Dye 3
25.	Acid Blue Dye 9
26.	N-Acetylprocainamide
27.	cis-Aconitic Acid
28.	Acrylamide
29.	Acrylic Acid
30.	Acrylonitrile
31.	Acyclovir
32.	Adenine
33.	Adenosine
34.	Adenosine 5'-diphosphate (ADP)
35.	Adenosine 5'-monophosphate (AMP)
36.	Adenosine 5'-monophosphate, Cyclic (cAMP)
37.	Adenosine 5'-triphosphate (ATP)
38.	Adipate

39.	Adipic Acid
40.	Adonitol
41.	Adrenaline
42.	Alanine
43.	Alkanesulfonates C2
44.	Alkanesulfonates C3
45.	Alkanesulfonates C4
46.	Alkanesulfonates C5
47.	Alkanesulfonates C6
48.	Alkanesulfonates C7
49.	Alkanesulfonates C8
50.	Alkanesulfonates C10
51.	Alkanesulfonates C12
52.	Alkanesulfonates C16
53.	Alkanesulfonates C18
54.	w-Alkoxy-3,5,7...-Polyoxaalkyl Sulfates (Alkyl Ether Sulfates)
55.	1-0-Alkyl-2-Acetyl-sn-Glycero-3-Phosphocholine (AGEPC)
56.	Alkylbenzenesulfonates
57.	Alkyl Sulfates
58.	Alkyl Sulfates C6
59.	Alkyl Sulfates C8
60.	Alkyl Sulfates C10
61.	Alkyl Sulfates C12
62.	Alkyl Sulfates C14
63.	Alkyl Sulfates C18
64.	Allantoin
65.	Allobarbitol
66.	Allylcyclopentenyl Barbaturic Acid
67.	Alphaphenal
68.	Alphaprodine
69.	Aluminum +3
70.	1-Aminoanthraquinone
71.	2-Aminoanthraquinone
72.	Aminoantipyrine
73.	p-Aminoazobenzene
74.	p-Aminobenzenesulfonic Acid
75.	3-Amino-4-Ethoxyacetanilide
76.	4-Amino-1-hydroxybutane-1,1-diphosphonic Acid Monosodium Salt Trihydrate
77.	1-Amino-2-Methylbenzene
78.	Aminomethylphosphonic Acid
79.	4-Amino-N ¹⁰ -methylpteroglutamic Acid
80.	1-Aminonaphthalene
81.	Aminooxyisobutyric Acid

82.	p-Aminophenol
83.	Aminopterin
84.	2-Aminopyrimidine
85.	Amitriptyline
86.	Ammonium
87.	Amobarbital
88.	AMP
89.	cAMP, (Cyclic AMP)
90.	AMPA
91.	Amphotericin
92.	Amphotericin B
93.	Amylocaine
94.	Androstene-3,17-dione
95.	Anhydroerythromycin A
96.	Anhydro-N-Demethyl-Erythromycin A
97.	Anhydro-Tetracycline
98.	Aniline
99.	o-Anisidine (Methoxyaniline)
100.	p-Anisidine
101.	Anisole
102.	Anthocyanidin
103.	Antimony (III)
104.	Antimony (V)
105.	Antipyrine
106.	Apoxytetracycline
107.	Apramycin
108.	Aprobarbital
109.	Arabinose
110.	Arabitol
111.	Arginine
112.	Arsenate
113.	Arsenic Acid
114.	Arsenite
115.	Arsenobetaine
116.	Arsenocholine
117.	Arsenous Acid
118.	Ascorbic Acid
119.	Asparagine
120.	Aspartame
121.	Aspartic Acid
122.	Aspirin
123.	Atropine
124.	5-Azacytidine
125.	Azelic Acid
126.	Azide
127.	Aztreonam (E)

128.	Aztreonam (Z)
129.	Aztreonam Desulfonated Open-ring Form
130.	Aztreonam Ethyl Ester (B)
131.	Aztreonam Open-ring Form (E)
132.	Aztreonam Open-ring Form (Z)
133.	Azure A (Asymmetrical Dimethyl-Thionin)
134.	Azure B
135.	Azure C
B	
136.	Barbital
137.	Barium
138.	Bentazon
139.	Benzaldehyde
140.	Benzamide
141.	Benzene
142.	Benzene Pentacarboxylic Acid
143.	Benzenesulfonic Acid
144.	Benzoate
145.	Benzocaine
146.	Benzoic Acid
147.	Benzonitrile
148.	Benzophenonetetracarboxylic Acid
149.	3,3',4,4'-Benzophenonetetracarboxylic Acid
150.	3,3',4,4'-Benzophenonetetracarboxylic Acid Dianhydride Ester
151.	3,3',4,4'-Benzophenonetetracarboxylic Acid Dimethyl Ester
152.	3,3',4,4'-Benzophenonetetracarboxylic Acid Monomethyl Ester
153.	3,3',4,4'-Benzophenonetetracarboxylic Acid Tetramethyl Ester
154.	3,3',4,4'-Benzophenonetetracarboxylic Acid Trimethyl Ester
155.	Benzoylgonine
156.	Benzphetamine
157.	Benzyl Alcohol
158.	Benzylamine
159.	Benzylammonium Ion
160.	Benzyl Carbamate
161.	Benzylcetyldimethylammonium chloride
162.	Benzyl Chloride
163.	Benzyltrimethylstearyl-ammonium chloride
164.	Benzyltrimethyltetradecyl-ammonium chloride

165.	Bethanechol
166.	Berberine
167.	BHA
168.	BHT
169.	Biliverdin
170.	Biotin
171.	Biphenyl
172.	Bis(acetylacetonato)cadmium(II)
173.	Bis(acetylacetonato)cobalt(II)
174.	Bis(acetylacetonato)copper(II)
175.	Bis(acetylacetonato)nickel(II)
176.	Bisdesethylchloroquine
177.	Borate
178.	Boronetrafluoride
179.	Bovine Cytochrome c
180.	Bovine Serum Albumin
181.	Bradykinin
182.	Bromate
183.	Bromite
184.	Bromobenzene
185.	Brucine
186.	Butabarbital
187.	Butacaine
188.	1,4-Butanediol
189.	1,4-Butanedione
190.	i-Butanol
191.	n-Butanol
192.	s-Butanol
193.	t-Butanol
194.	Butanone
195.	Butethal
196.	Butirosin A
197.	Butirosine
198.	Butylamine
199.	Butylated Hydroxyanisole
200.	Butylated Hydroxytoluene
201.	Butyric Acid
202.	Butyrophenone
C	
203.	C14-AGEPC
204.	C16-AGEPC
205.	C18-AGEPC
206.	Cadmium
207.	Cadmium Phenanthroline
208.	Caffeic Acid

209.	Caffeine
210.	Calciferol
211.	Calcium
212.	Calmodulin
213.	Canine Cytochrome c
214.	Cannabinol
215.	Capric Acid
216.	Caproic Acid
217.	Caprylic Acid
218.	Carbamazepine
219.	Carbonate
220.	2-Carboxy-D-Arabinitol-1,5-Bisphosphate
221.	2-Carboxy-D-Arabinitol-1-Phosphate
222.	2-Carboxy-D-Arabinitol-5-Phosphate
223.	2-Carboxy-D-Ribitol-1,5-Bisphosphate
224.	Carvacrole
225.	Carvacrol (5-Isopropyl-2-Methyl-Phenol)
226.	Catechin
227.	CBZ-1-Azetidinone
228.	CBZ-1-Azetidinone-TBA Salt
229.	CBZ-1-Serinamide
230.	CBZ-1-Serinamide Mesylate
231.	CBZ-1-Serine
232.	CBZ-1-Threoninamide
233.	CBZ-1-Threonine
234.	Cesium
235.	Chenodeoxycholate, Glycine and Taurine-conjugated
236.	Chenodeoxycholic Acid
237.	Chimeric Ribozyme 32mer
238.	Chimeric Ribozyme 38mer
239.	Chloramphenicol
240.	Chloranil (Tetrachloro-1,4-Benzodionin)
241.	Chlorate
242.	Chlordiazepoxide
243.	Chlordimeform
244.	Chlorhexidine Gluconate
245.	Chlorite
246.	Chloroacetic Acid
247.	p-Chloroaniline
248.	Chlorobenzene
249.	p-Chlorobenzenesulfonate
250.	3-Chlorobenzohydroxamic Acid
251.	3-Chlorobenzoic Acid

252.	2-Chloroethylethyl-Methylsulfonium Ion
253.	Chloroform
254.	Chlorogenic Acid
255.	4-Chloro-2-Methylphenoxy-Acetic Acid
256.	4-(4-Chloro-2-methylphenoxy)-butanoic Acid
257.	2-(4-Chloro-2-methylphenoxy)-Propionic Acid
258.	p-Chloronitrobenzene
259.	2-Chlorophenol
260.	3-Chlorophenol
261.	4-Chlorophenol
262.	2-Chlorophenoxyacetic Acid
263.	4-Chlorophenoxyacetic Acid
264.	2-Chloropropanoic Acid
265.	3-Chloropropanoic Acid
266.	Chloroquine
267.	Chlor-l-succinate
268.	Chlor-3-succinate
269.	7-Chlorotetracycline
270.	8-Chlorotheophylline
271.	Chloroxylene
272.	Chlorphenesin Carbamate
273.	Chlorpheniramine
274.	Chlorpromazine
275.	Chlorpropamine
276.	Chlortetracycline
277.	4-epi-Chlortetracycline
278.	keto-4-Chlortetracycline
279.	keto-4-epi-Chlortetracycline
280.	Cholate, Glycine and Taurine-conjugated
281.	Cholcalciferol
282.	Cholic Acid
283.	Choline
284.	Chromate
285.	Chromium
286.	Chromosomal DNA
287.	Cinchonidine
288.	Cinchonine
289.	Ciprofloxacin
290.	Cisplatin
291.	Citrate
292.	Citric Acid
293.	Clomesone
294.	Clonazepam
295.	Clonidine

296.	CMP
297.	cCMP, (Cyclic CMP)
298.	Cobalt
299.	Cocaine
300.	Codeine
301.	Coenzyme Q10
302.	Colchicine
303.	Conalbumin
304.	Concanavalin A
305.	Contrathion
306.	Copper
307.	Copper EDTA
308.	Cortisone
309.	o-Coumaric Acid
310.	p-Coumaric Acid
311.	Coumarin
312.	Crayfish Carboxypeptidase B Tryptic Digest
313.	Creatinine
314.	m-Cresol
315.	o-Cresol
316.	p-Cresol
317.	Cyanate
318.	Cyanide
319.	Cyanocobalamin
320.	Cyclamic Acid
321.	Cyclic P3
322.	Cyclodisone
323.	1,2-Cyclohexanedione
324.	1,2-Cyclohexanedione Dioxime
325.	Cyclohexanedione Monoxime
326.	Cyclohexanone Monoxime
327.	Cyproheptadine
328.	Cysteine
329.	Cytidine
330.	Cytidine 5'-diphosphate (CDP)
331.	Cytidine 5'-monophosphate (CMP)
332.	Cytidine 5'-monophosphate, Cyclic (cCMP)
333.	Cytidine 5'-triphosphate (CTP)
334.	Cytochrome
335.	Cytosine
D	
336.	2,4-D
337.	Daidzein
338.	Dansyl-Alanine

339.	Dansyl-Arginine
340.	Dansyl-Asparagine
341.	Dansyl-Glutamine
342.	Dansyl-Glycine
343.	Dansyl-Leucine
344.	Dansyl-Lysine
345.	Dansyl-Methionine
346.	Dansyl-Phenylalanine
347.	Dansyl-Serine
348.	Dansyl-Tryptophan
349.	Dansyl-Tyrosine
350.	Danthron
351.	Dapsone
352.	DCTA (1,2-Diaminocyclohexane-Tetraacetic Acid)
353.	Decane
354.	1-Decanol
355.	2-Decanone
356.	Decarbamoylgonyautoxin I
357.	Decarbamoylgonyautoxin II
358.	Decarbamoylgonyautoxin III
359.	Decarbamoylgonyautoxin IV
360.	Decarbamoylneosaxitoxin
361.	Decarbamoylsaxitoxin
362.	Dehydroascorbic Acid
363.	Dehydroepiandrosterone
364.	Dehydroisoascorbic Acid
365.	Demeclocycline
366.	Demethylchlortetracycline
367.	6-Demethylchlortetracycline
368.	N-Demethylerythromycin A
369.	N-Demethylerythromycin A Enol Ether
370.	Deoxycholate, Glycine and Taurine-conjugated
371.	Deoxycholic Acid
372.	2-Deoxyribose
373.	Deprotected AAF-Modified Oligonucleotide 6mer
374.	Deprotected AAF-Modified Oligonucleotide 9mer
375.	Deprotected AAF-Modified Oligonucleotide 12mer
376.	Deprotected Oligonucleotide 4mer
377.	Deprotected Oligonucleotide 6mer
378.	Deprotected Oligonucleotide 8mer
379.	Deprotected Oligonucleotide 9mer
380.	Deprotected Oligonucleotide 10mer

381.	Deprotected Oligonucleotide 11mer
382.	Deprotected Oligonucleotide 12mer
383.	Deprotected Oligonucleotide 15mer
384.	Deprotected Oligonucleotide 16mer
385.	Deprotected Oligonucleotide 17mer
386.	Deprotected Oligonucleotide 20mer
387.	Deprotected Oligonucleotide 21mer
388.	Deprotected Oligonucleotide 22mer
389.	Deprotected Oligonucleotide 24mer
390.	Deprotected Oligonucleotide 28mer
391.	Deprotected Oligonucleotide 32mer
392.	Deprotected Oligonucleotide 36mer
393.	Deprotected Oligonucleotide 44mer
394.	Deprotected Oligonucleotide 50mer
395.	Desethylchloroquine
396.	Desethylhydroxychloroquine
397.	N-DesmethyImphenytoin
398.	N-DesmethyImethsuximide
399.	DespropanylFentanyl
400.	Desulfonated Aztreonam
401.	Detritylated Oligonucleotide 8mer
402.	Detritylated Oligonucleotide 16mer
403.	Detritylated Oligonucleotide 20mer
404.	Dexamethasone
405.	1,1'-Diacetylferrocene
406.	Diallylbarbituric Acid
407.	Diallylmethylsulfonium
408.	1,2-Diaminocyclohexane-Tetraacetic Acid
409.	cis-Diamminedichloroplatinum (II)
410.	Dianisidine (Dimethoxybenzidine)
411.	o-Dianisidine
412.	Diazepam
413.	Dibucane
414.	Dibutylamine
415.	Dibutylphosphate
416.	Dicamba
417.	1,2-Dichlorobenzene
418.	1,4-Dichlorobenzene
419.	2,5-Dichlorobenzene Sulfonic Acid
420.	3,3'-Dichlorobenzidine
421.	2,4-Dichlorobenzohydroxamic Acid
422.	3,6-Dichloro-2-Methoxybenzoic Acid
423.	2,3-Dichlorophenol
424.	2,4-Dichlorophenol
425.	2,5-Dichlorophenol

426.	2,6-Dichlorophenol
427.	3,4-Dichlorophenol
428.	3,5-Dichlorophenol
429.	2,4-Dichlorophenoxyacetic Acid (2,4-D)
430.	4-(2,4-Dichlorophenoxy)butanoic Acid
431.	2-(2,4-Dichlorophenoxy)propionic Acid
432.	2,2-Dichloropropanoic Acid
433.	2,3-Dichloropropanoic Acid
434.	Dicyanodiamide
435.	Didrate
436.	Diethanolamine
437.	Diethylamine
438.	Diethylene Glycol
439.	Diethylenetriaminepentaacetic Acid
440.	Diethylmethylsulfonium
441.	Diethylphthalate
442.	Diethylstilbestrol
443.	Diethyltryptamine
444.	Difacinone
445.	Diffluoro(Acetylacetonato) Boron(III)
446.	Diglycolate
447.	6,7-Dihydro-7-Hydroxy-1-Hydroxymethyl-5h-Pyrrolizine
448.	Dihydroxyacetone
449.	1,2 Dihydroxybenzene
450.	Dilaudid
451.	Dilevalol
452.	Diltiazem
453.	Dimethoxybenzidine (Dianisidine)
454.	Dimethylallylsulfonium Ion
455.	Dimethylamine
456.	Dimethylarsinic Acid
457.	Dimethylarsinous Acid
458.	1-3-Dimethylbarbituric Acid
459.	2,5-Dimethylbenzenesulfonic Acid
460.	1,1-Dimethylguanidine
461.	Dimethyl-2-Hydroxyethyl-Sulfonium
462.	2,3-Dimethyl-1-(4-Methyl-Phenyl)-3-Pyrazolin-5-one
463.	Dimethylphthalate
464.	Dimethylpropylsulfonium Ion
465.	Dimethylselenonium Iodide
466.	Dimethylthiourea
467.	2,4-Dinitrobenzenesulfonic Acid
468.	3,5-Dinitro-o-cresol
469.	2,4-Dinitrophenol
470.	2,6-Dinitrophenol

471.	2,4-Dinitrotoluene-3-sulfonic Acid
472.	Dipeptides
473.	Diphenoxylate
474.	Diphenylguanidine
475.	Diphenylhydantoin
476.	Diphenylhydramine
477.	Dipropyltryptamine
478.	Diquat
479.	Disodium Edetate
480.	2,6-Di-tert-butyl-p-cresol
481.	Dithionate
482.	Dithiothreitol
483.	DMT Protected AAF Modifier Oligonucleotide 12mer
484.	DMT Protected Oligonucleotide 4mer
485.	DMT Protected Oligonucleotide 6mer
486.	DMT Protected Oligonucleotide 8mer
487.	DMT Protected Oligonucleotide 9mer
488.	DMT Protected Oligonucleotide 10mer
489.	DMT Protected Oligonucleotide 12mer
490.	DMT Protected Oligonucleotide 15mer
491.	DMT Protected Oligonucleotide 16mer
492.	DMT Protected Oligonucleotide 17mer
493.	DMT Protected Oligonucleotide 20mer
494.	DMT Protected Oligonucleotide 21mer
495.	DMT Protected Oligonucleotide 24mer
496.	DMT Protected Oligonucleotide 25mer
497.	DMT Protected Oligonucleotide 28mer
498.	DMT Protected Oligonucleotide 30mer
499.	DMT Protected Oligonucleotide 36mer
500.	DMT Protected Oligonucleotide 44mer
501.	DMT Protected Oligonucleotide 50mer
502.	DNA Fragments of p-BR322
503.	DNA Oligonucleotide 10mer
504.	DNA Oligonucleotide 11mer
505.	DNA Oligonucleotide 12mer
506.	DNA Oligonucleotide 16mer
507.	DNA Oligonucleotide 20mer
508.	DNA Oligonucleotide 21mer
509.	DNA Oligonucleotide 32mer
510.	DNP-Alanine
511.	DNP-Arginine
512.	DNP-Asparagine
513.	DNP-Glutamic Acid
514.	DNP-Glutamine
515.	DNP-Leucine

516.	DNP-Lysine
517.	DNP-Methionine
518.	DNP-Phenylalanine
519.	DNP-Proline
520.	DNP-Serine
521.	DNP-Threonine
522.	DNP-Tryptophan
523.	DNP-Tyrosine
524.	Dodecane
525.	Dodecanesulfonate
526.	Dodecyl Sulfate
527.	Dog Albumin
528.	l-Dopa
529.	Dopamine
530.	Dowicil-75
531.	Doxapram
532.	Doxycycline
533.	DP2 Oligosaccharide
534.	DP5 Oligosaccharide
535.	DP10 Oligosaccharide
536.	DP15 Oligosaccharide
537.	Dyphylline

E

538.	Edrophonium
539.	EDTA (Ethylenediaminetetraacetic Acid)
540.	EDTA, Copper
541.	EDTA, Iron
542.	Ellagic Acid
543.	Enalapril Maleate
544.	Endothall Acid
545.	Enoxacin
546.	Enthanic Acid
547.	Ephedrine
548.	Epianhydrotetracycline
549.	Epicatechin
550.	4-Epidoxycycline
551.	6-Epidoxycycline
552.	4,6-Epidoxycycline
553.	Epitetracycline
554.	Equine Cytochrome c
555.	Ergocristine
556.	Ergocristinine
557.	α -Ergocryptine
558.	α -Ergocryptinine
559.	Ergonovine

560.	Ergonovine
561.	Ergotamine
562.	Ergotaminine
563.	Erioglaucine
564.	Erythralosamine
565.	Erythritol
566.	Erythromycin
567.	Erythromycin A
568.	Erythromycin A Enol Ether
569.	Erythromycin B
570.	Erythromycin C
571.	Erythromycin Related Compound N
572.	Estrone
573.	Ethanol
574.	Ethanolamine
575.	Ethonitazene
576.	Ethosuximide
577.	Ethoxylated Alcohols C6 to C18
578.	Ethylamine
579.	9-Ethyl-3-Aminocarbazole
580.	Ethylbenzene
581.	Ethyl Benzoate
582.	<i>O</i> -Ethyl <i>S</i> -[2[(diisopropylamino)ethyl]isopropylphosphonothiolate
583.	<i>O</i> -Ethyl <i>S</i> -[2[(diisopropylamino)ethyl]propylphosphonothiolate
584.	<i>O</i> -Ethyl <i>S</i> -[2[(dimethylamino)ethyl]- <i>n</i> -propylphosphonothiolate
585.	<i>O</i> -Ethyl <i>S</i> -[2[(dimethylamino)ethyl]propylphosphonothiolate
586.	Ethylene Cyanohydrin
587.	Ethylenediamine
588.	Ethylenediaminetetraacetic Acid
589.	Ethylene Glycol
590.	1-Ethylguanidine
591.	2-Ethylhexylsalicylate
592.	2-Ethylhexyl <i>trans</i> -4-methoxycinnamate
593.	Ethyl Homocholine
594.	Ethylmorphine
595.	<i>N</i> -Ethylornicotine
596.	5-Ethyl-5- <i>p</i> -Tolylbarbituric Acid
597.	Etonitazene
598.	Eugenol
599.	Everninomycin, SCH 27899

F

600.	FAD
------	-----

601.	Fenfluramine
602.	Fentanyl
603.	Fentanyl Homologues and Analogues
604.	Ferrocene
605.	Ferulic Acid
606.	Flavin Adenine Dinucleotide
607.	Fluorescein Labeled Oligonucleotide 12mer
608.	Fluoride
609.	Fluoroacetic Acid
610.	Fluoroborate
611.	5-Fluorouracil
612.	Flurazepam
613.	Folic Acid
614.	Formaldehyde
615.	Formate
616.	Formic Acid
617.	Fortimicin A
618.	Fortimicin B
619.	Fructose
620.	Fructose-1,6-Diphosphate
621.	Fructose-6-Phosphate
622.	Fucose
623.	Fumaric Acid
624.	Furfural

G

625.	Galactitol
626.	Galactosamine
627.	Galactose
628.	Galacturonic Acid
629.	Gallic Acid
630.	Genistein
631.	Gentamicin C1
632.	Gentamicin C1a
633.	Gentamicin C2
634.	Gentamicin C2a
635.	Gentisic Acid
636.	Gerall-115
637.	Gibberellic Acid
638.	Gibberellin A ₇
639.	Glucaric Acid
640.	Glucobrassicin
641.	Gluconapin
642.	Gluconasurtiin
643.	Gluconic Acid
644.	Glucosamine

645.	Glucosamine-1-phosphate
646.	Glucosamine-6-phosphate
647.	Glucose-6-Phosphate
648.	Glucotropaeolin
649.	Glucuronic Acid
650.	Glufosinate Ammonium
651.	Glutamate
652.	Glutamic Acid
653.	Glutamine
654.	g-Glutamylglutamine
655.	Glutarate
656.	Glutaric Acid
657.	Glutathione
658.	Glutethimide
659.	Glycerine
660.	Glycerol
661.	Glycine
662.	Glycine-conjugated, Chenodeoxycholate
663.	Glycine-conjugated Cholate
664.	Glycine-conjugated, Deoxycholate
665.	Glycine-conjugated, Lithocholate
666.	Glycine-conjugated, Ursodeoxycholate
667.	Glycolic Acid
668.	Glycylmethoxyacetic Acid
669.	Glycylsarcosine
670.	Glycyrrhizic Acid
671.	Glydant
672.	Glyoxal
673.	Glyoxylate
674.	Glyphosate
675.	GMP
676.	cGMP, (Cyclic GMP)
677.	Gonyautoxin I
678.	Gonyautoxin II
679.	Gonyautoxin III
680.	Gonyautoxin IV
681.	Gonyautoxin VIII
682.	Guaiacol
683.	Guanidine
684.	Guanine
685.	Guanosine
686.	Guanosine 5'-diphosphate (GDP)
687.	Guanosine 5'-monophosphate (GMP)
688.	Guanosine 5'-monophosphate, Cyclic (cGMP)
689.	Guanosine 5'-triphosphate (GTP)

H	
690.	HEDTA (Hydroxyethylene-Diaminetriacetate)
691.	Heptane
692.	1-Heptanol
693.	2-Heptanone
694.	Heptanophenone
695.	Heptapolyphosphate
696.	Hexabarbital
697.	Hexamethylenetetramine
698.	1-Hexanol
699.	2-Hexanone
700.	Hexanophenone
701.	Hexapolyphosphate
702.	Hexylamine
703.	Histamine
704.	Histidine
705.	Hoelan
706.	Hyaluronic Acid Decasaccharide
707.	Hyaluronic Acid Dodecasaccharide
708.	Hyaluronic Acid Hexasaccharide
709.	Hyaluronic Acid Octasaccharide
710.	Hyaluronic Acid Tetrasaccharide
711.	Hydrastine
712.	Hydrazine
713.	Hydrocortisone
714.	Hydrogen Sulfide
715.	m-Hydroxybenzaldehyde
716.	Hydroxybenzoate
717.	m-Hydroxybenzoic Acid
718.	o-Hydroxybenzoic Acid
719.	p-Hydroxybenzoic Acid
720.	a-Hydroxybutyric Acid
721.	g-Hydroxybutyric Acid
722.	Hydroxychloroquine
723.	9-((2-Hydroxyethoxy)methyl)Guanine
724.	2-Hydroxyethylethylmethylsulfonium
725.	Hydroxyethyltheophylline
726.	4-Hydroxy-Glucobrassicin
727.	5-Hydroxyisoquinol
728.	7-Hydroxymethotrexate
729.	2-Hydroxy-4-methoxybenzophenone
730.	Hydroxymethyl-furfural
731.	5-Hydroxy-2-Pyrrolidinone
732.	2-Hydroxyquinoline
733.	4-Hydroxyquinoline

734.	8-Hydroxyquinoline
735.	Hyoscyamine
736.	Hypophosphite
737.	Hypoxanthine
I	
738.	Ibuprofen
739.	Imipramine
740.	Imminodiacetic Acid
741.	Indol-3-Carboxaldehyde
742.	Indole-3-Acetic Acid
743.	Indomethacin
744.	Inosine
745.	Inositol
746.	Inositol Hexaphosphate
747.	Inositol Pentaphosphate
748.	Inositol Phosphate
749.	Inositol Tetrphosphate
750.	Inositol Triphosphate
751.	Insulin
752.	Iodate
753.	Iodide
754.	Iron
755.	Iron Diethyletriaminepentaacetic Acid
756.	Iron Ethylenediaminetetraacetic Acid (EDTA)
757.	Iron Nitritotriacetic Acid
758.	Iron Phenanthroline
759.	Isoascorbic Acid
760.	3-Isobutyl-1-Methylxanthine
761.	R,S-2-(4-Isobutylphenyl)-Propionic Acid
762.	Isobutyric Acid
763.	Isochlortetracycline
764.	Isocitric Acid
765.	Isoleucine
766.	Isomalt
767.	Isomaltose
768.	Isomaltotriose
769.	Isopropanol
770.	3-Isopropyl-(1 <i>H</i>)-2,1,3-Benzo-thiadiazin-4(3 <i>H</i>)-one 2,2 Dioxide
771.	Isopropyl-6-Methylphenol
772.	2-Isopropyl-5-Methylphenol (Thymol)
773.	5-Isopropyl-2-Methylphenol (Carvacrol)
774.	Isopropylmethyl Phosphonic Acid
775.	Isoquinoline-N-Oxide
776.	Isovaleric Acid

J, K

777.	Jacobine
778.	Jacozine
779.	Kaempferol
780.	Kanamycin
781.	Kanamycin A
782.	α -Ketobutyric Acid
783.	α -ketocaproic Acid
784.	α -Ketoglutaric Acid
785.	α -Ketoisocaproic Acid
786.	α -Ketoisovaleric Acid
787.	Ketomalonic
788.	6-Keto-Prostaglandin F _{1a}

L

789.	Labetalol
790.	Lactate
791.	Lactate Dehydrogenase
792.	Lactic Acid
793.	Lactose
794.	Lauric Acid
795.	Lead
796.	Leucine
797.	Levorphanol
798.	LHC II
799.	Lidocaine
800.	Linamarin
801.	Linoleic Acid
802.	Linolenic Acid
803.	Liothyronine (T3)
804.	Lithium
805.	Litocholate, Glycine and Taurine-conjugated
806.	Litocholic Acid
807.	Lividomycin A
808.	Lomefloxacin
809.	Lutidine
810.	Lymecycline
811.	Lysine
812.	Lysozyme
813.	Lysozyme Tryptic Digest Reduced and s-Carboxymethylated

M

814.	Magnesium
815.	Malate
816.	Maleic Acid

817.	Malic Acid
818.	d-Malic Acid
819.	l-Malic Acid
820.	Malonate
821.	Malonic Acid,
822.	Maltodecaose
823.	Maltoheptaose
824.	Maltohexaose
825.	Maltononaose
826.	Maltooctaose
827.	Maltopentaose
828.	Maltose
829.	Maltotetraose
830.	Maltotriose
831.	Manganese
832.	Manganese (III) Protoporphyrin
833.	Mannitol
834.	Mannose
835.	Meclizine
836.	Mefenamic Acid
837.	Mepenzolate
838.	Meperidine
839.	Mephensin
840.	Mephenytoin
841.	Mephobarbital
842.	Mepivacine
843.	Mercaptoacetic Acid
844.	Mescaline
845.	Metacycline
846.	Metadrenaline
847.	Metaphosphate
848.	Methacycline
849.	Methadone
850.	Methagalone
851.	Metham Sodium
852.	Methanol
853.	Methionine
854.	Methocarbamol
855.	Methotrexate
856.	Methoxamine
857.	Methoxyacetic Acid
858.	Methoxyaniline (o-Anisidine)
859.	4-Methoxybenzohydroxamic Acid
860.	2-Methoxyethanol
861.	2-Methoxyphenol
862.	Methsuxamide

863.	Methylamine
864.	N-Methylaniline
865.	Methylarsonic Acid
866.	Methylarsonous Acid
867.	p-Methylbenzenesulfonic Acid
868.	Methyl Benzoate
869.	4-Methylbenzohydroxamic Acid
870.	B-Methylcholine
871.	1-Methyl-2,4-Diaminobenzene
872.	Methyl (+)2-(4-(2,4dichlorophenoxy)) propionate (diclofopmethyl)
873.	Methyldopa
874.	Methylene Blue
875.	4,4'-Methylenedianiline
876.	4,4'-Methylenedianiline Bisnadimide
877.	4,4'-Methylenedianiline Mononadamide
878.	Methyl-furfural
879.	N-Methylglucamine
880.	1-Methylguanidine
881.	Methyl Isocyanate
882.	Methylmalonic Acid
883.	N-Methylmorpholine Oxide
884.	N-Methylphenazone
885.	Methylphenidate
886.	4-Methyl-m-Phenylenediamine
887.	3-Methylxanthine
888.	Minocycline
889.	Molybdate
890.	Monobutyl Phosphate
891.	Monochloroacetic Acid
892.	Monofluorophosphate
893.	Monohydrochloride
894.	Monomethylarsonic Acid
895.	Monosodium Glutamate
896.	Morphine
897.	Morpholine
898.	Muconic Acid
899.	Myoglobin Myoglobin Tryptic Digest
900.	Myristic Acid

N

901.	Nadolol
902.	Nalidixic Acid
903.	Nalorphine
904.	Naloxone

905.	6-B-Naltrexol
906.	Naphthalene
907.	2-[(1-Naphthalenylamino)-Carbonyl] Benzoic Acid
908.	1-Naphthylamine
909.	Naproxen
910.	Naptalam
911.	n-Butanol
912.	Neomycin B
913.	Neosaxitoxin
914.	Neostigmine
915.	Niacin
916.	Niacinamide
917.	Nickel
918.	Nicotine
919.	Nicotinic Acid
920.	Nikethamide
921.	Nitrate
922.	Nitrilotriacetic Acid
923.	Nitrite
924.	p-Nitroaniline
925.	Nitrobenzene
926.	m-Nitrobenzenesulfonate
927.	p-Nitrobenzenesulfonic Acid
928.	3-Nitro-4-hydroxyphenylarsonic Acid
929.	Nitromethane
930.	p-Nitrophenol
931.	Nonane
932.	1-Nonanol
933.	2-Nonanone
934.	Noradrenaline
935.	5-Norbornene-2,3-Dicarboxylic Acid Anhydride
936.	5-Norbornene-2,3-Dicarboxylic Acid Dimethyl Ester
937.	5-Norbornene-2,3-Dicarboxylic Acid Monomethyl Ester
938.	Norfloxacin
939.	Norleucine
940.	Normetadrenaline
941.	N-Normorphine

O

942.	Octadecylsulfate
943.	Octane
944.	1-Octanol
945.	2-Octanone

946.	Ofloxacin
947.	Oleandomycin
948.	Oleic Acid
949.	Oligoribonucleotides
950.	Ornithine
951.	Orthophosphate
952.	Ovalbumin
953.	Oxalacetic Acid
954.	Oxalate
955.	Oxalic Acid
956.	Oxaloacetic Acid
957.	Oxalosuccinic Acid
958.	Oxazepam
959.	Oxobis(Acetylacetonato) Vanadium(IV)
960.	Oxolinic Acid
961.	Oxytetracycline

P

962.	Palmitic Acid
963.	Panthenol
964.	Papaverine
965.	Paraquat
966.	Paraxanthine
967.	Patent Blue VF
968.	PCP N-Ethyl Analog
969.	PEG 600
970.	Pelargonic Acid
971.	Pemoline
972.	2,3,4,5,6-Pentachlorophenol
973.	Pentaerythritol
974.	2-Pentanone
975.	Pentapolyphosphate
976.	Pentazocine
977.	Pentobarbital
978.	Pentylamine
979.	Peptides
980.	Perchlorate
981.	Permanganate
982.	Phenacetate
983.	Phenacetin
984.	Phencyclidine
985.	Phenelzine
986.	Phenobarbital
987.	Phenol
988.	2-Phenylaminonapthalene

989.	4-Phenylbutanoate
990.	Phentermine
991.	Phenylacetate
992.	Phenylalanine
993.	Phenylbutylammonium Ion
994.	10-Phenyldecanohydroxamic Acid
995.	m-Phenylenediamine
996.	o-Phenylenediamine
997.	p-Phenylenediamine
998.	2-Phenylethanol
999.	Phenylethylammonium Ion
1000.	Phenylethylmalonamide
1001.	Phenylglucuronide
1002.	6-Phenylhexanohydroxamic Acid
1003.	Phenylmercapturate
1004.	9-Phenylnonanohydroxamic Acid
1005.	9-Phenylnonanoic Acid
1006.	8-Phenyloctanohydroxamic Acid
1007.	5-Phenylpentanoate
1008.	3-Phenylpropanoate
1009.	Phenylpropylammonium Ion
1010.	Phenylsulfate
1011.	Phorwite RKH
1012.	Phosphate
1013.	Phosphate Oligomer P1
1014.	Phosphate Oligomer P2
1015.	Phosphate Oligomer P3
1016.	Phosphate Oligomer P4
1017.	Phosphate Oligomer P5
1018.	Phosphate Oligomer P6
1019.	Phosphate Oligomer P7
1020.	Phosphate Oligomer P8
1021.	Phosphate Oligomer P9
1022.	Phosphate Oligomer P10
1023.	Phosphate Oligomer P11
1024.	Phosphate Oligomer P12
1025.	Phosphate Oligomer P13
1026.	Phosphatidylcholine
1027.	Phosphatidylethanolamine
1028.	Phosphite
1029.	Phosphoarginine
1030.	3-Phosphoglyceric Acid
1031.	1-Phosphohistidine
1032.	3-Phosphohistidine
1033.	N-(Phosphonomethyl)glycine
1034.	N-Phosphonomethyliminodiacetic Acid
1035.	Phosphorothioate 28mer

1036.	Phosphorothioate deoxycytidine homopolymer 28mer
1037.	Phosphoserine
1038.	Phosphothreonine
1039.	Phosphotyrosine
1040.	Photosystem II Light-Harvesting Complex (LHC II)
1041.	Phthalic Acid
1042.	Phylloquinone
1043.	Physostigmine
1044.	Phytic Acid
1045.	Picramic Acid
1046.	Pimelic Acid
1047.	Pipemidic Acid
1048.	Pipenzolate
1049.	Piromidic Acid
1050.	Plasmid DNA
1051.	Polyethylene Glycol
1052.	Polyethylene Glycol 600
1053.	Polysulfide
1054.	Polythiazide
1055.	Potassium
1056.	Potassium Tetrathionate
1057.	Prednisolone
1058.	Prednisone
1059.	Pregnenolone
1060.	Pregnenolone-3-sulfate
1061.	Primidone
1062.	Procaine
1063.	Progesterone
1064.	Progoitrin
1065.	Proline
1066.	L-Proline 3,4- ³ H(N),histidyl-3- ³ H(N)
1067.	Promazine
1068.	Promethazine
1069.	Propanoic Acid
1070.	1-Propanol
1071.	2-Propanol
1072.	2-Propanone
1073.	Propiomazine
1074.	Propionate
1075.	Propionic Acid
1076.	Propiophenone
1077.	Propranolol
1078.	Propylamine
1079.	Propylene Glycol
1080.	Prostacyclin (PGI ₂)

1081.	Protected DNA Oligomers
1082.	Proteins
1083.	Protocatechuic Acid
1084.	Protoporphyrin IX
1085.	PTH-Alanine
1086.	PTH-Arginine
1087.	PTH-Asparagine
1088.	PTH-Glutamic Acid
1089.	PTH-Glycine
1090.	PTH-Leucine
1091.	PTH-Lysine
1092.	PTH-Phenylalanine
1093.	PTH-Tryptophan
1094.	PTH-Tyrosine
1095.	PTH-Valine
1096.	Pyridine
1097.	Pyridinoline
1098.	Pyridostigmine
1099.	Pyridoxine
1100.	Pyrilamine
1101.	L-Pyroglutamyl-L-histidyl-L-proline amide
1102.	Pyrophosphate
1103.	2-Pyrrolidinone
1104.	Pyruvate
1105.	Pyruvate Dehydrogenase
1106.	Pyruvic Acid

Q, R

1107.	Quercetin
1108.	Quinoline-N-Oxide
1109.	Quinolones
1110.	Ranitidine
1111.	RB
1112.	Reserpine
1113.	Resveratrol
1114.	Retinol
1115.	Retinol Acetate
1116.	Retrorsine
1117.	Retrorsine N-Oxide
1118.	Rhamnose
1119.	Ribitol
1120.	Riboflavin
1121.	Riboflavin 5'-Phosphate
1122.	Ribonuclease A
1123.	Ribose

1124.	Ribose-5-Phosphate
1125.	Ribostamycin
1126.	Ribozyme 32mer, Chimeric
1127.	Ribozyme 38mer, Chimeric
1128.	Ribulose-1,5-Diphosphate
1129.	RNA
1130.	Rolitetraacycline
1131.	Rubidium
1132.	Ruthenium Phenanthroline

S

1133.	Saccharin
1134.	Salicylamide
1135.	Salicylic Acid
1136.	Saligenin
1137.	Saxitoxin
1138.	Scopolamine
1139.	Scopoletin
1140.	Secobarbital
1141.	Sedoheptulose-1,7-Diphosphate
1142.	Selenate
1143.	Selenic Acid
1144.	Selenite
1145.	Selenocystine
1146.	Selenoethionine
1147.	Selenohomocystine
1148.	Selenomethionine
1149.	Selenous Acid
1150.	Senecionine
1151.	Senecionine N-Oxide
1152.	Seneciphylline
1153.	Serine
1154.	Shikimic Acid
1155.	Sialic Acid, (N-Acetyl Neuraminic Acid)
1156.	Silicate
1157.	Silver Iodide
1158.	Sinalbin
1159.	Sinapic Acid
1160.	Sinigrin
1161.	Sisomicin
1162.	SKF 78729A
1163.	Slaframine
1164.	Sodium
1165.	Sorbic Acid
1166.	Sorbitol

1167. Soybean Trypsin Inhibitor
1168. Stearic Acid
1169. Streptomycin
1170. Strontium
1171. Strychnine
1172. Suberic Acid
1173. Succinate
1174. Succinic Acid
1175. Succinimide
1176. Succinylsulfathiazole
1177. Sucrose
1178. Sulfacetamide
1179. Sulfadiazine
1180. Sulfadimethoxine
1181. Sulfaguanidine
1182. Sulfamate
1183. Sulfamerazine
1184. Sulfamethazine
1185. Sulfamethizole
1186. Sulfamethoxazole
1187. Sulfanilamide
1188. Sulfanilic Acid
1189. Sulfanilylbenzamide
1190. Sulfapyridine
1191. Sulfate
1192. Sulfathiazole
1193. Sulfide
1194. Sulfisoxazole
1195. Sulfite
1196. Sulfite as Alphahydroxymethylsulfonate
1197. <i>N</i> -Sulfocarbamoyl Toxin B1
1198. <i>N</i> -Sulfocarbamoyl Toxin B2
1199. <i>N</i> -Sulfocarbamoyl Toxin C1
1200. <i>N</i> -Sulfocarbamoyl Toxin C2
1201. <i>N</i> -Sulfocarbamoyl Toxin C3
1202. <i>N</i> -Sulfocarbamoyl Toxin C4
1203. Symmetrical Dimethylthionin
1204. Syringaldehyde
1205. Syringic Acid
T
1206. T3
1207. T4
1208. Tartaric Acid
1209. <i>d</i> -Tartaric Acid
1210. <i>l</i> -Tartaric Acid

1211. Tartaric Acid
1212. Taurine
1213. Taurine-conjugated Chenodeoxycholate
1214. Taurine-conjugated Cholate
1215. Taurine-conjugated Deoxycholate
1216. Taurine-conjugated Lithocholate
1217. Taurine-conjugated Ursodeoxycholate
1218. ^{99m} Tc Diethylenetriaminepentaacetic Acid
1219. ^{99m} Tc Glucoheptonate
1220. ^{99m} Tc Methylene Diphosphonate
1221. ^{99m} Tc Pertechnetate
1222. <i>o</i> -Terphenyl
1223. Testosterone
1224. Tetrabutylammonium Chloride
1225. Tetracaine
1226. 2,3,5,6-Tetracarboxybenzene Sulfonic Acid
1227. Tetrachloro-1,4-Benzodion (Chloranil)
1228. 2,3,4,5-Tetrachlorophenol
1229. 2,3,4,6-Tetrachlorophenol
1230. 2,3,5,6-Tetrachlorophenol
1231. Tetracycline
1232. Tetradecane
1233. Tetradecyl Sulfate
1234. Tetraethylammonium Chloride
1235. Tetramethylammonium Hydroxide
1236. Tetramethylarsonium Iodide
1237. Tetramethylarsonium Ion
1238. Tetramethylenemethylsulfonium Ion
1239. 1,1,3,3-Tetramethylguanidine
1240. Tetrapentylammonium Chloride
1241. Tetrapolyphosphate
1242. Tetrapropylammonium Chloride
1243. Tetrathionate
1244. Theaflavin
1245. Theaflavin-3'-digallate
1246. Theaflavin-3-gallate
1247. Theaflavin-3'-gallate
1248. Theafulvins
1249. Thebaine
1250. Theophylline
1251. Thiabendazole
1252. Thiamine

1253. Thiamine Diphosphate
1254. Thiamine Monophosphate
1255. Thiamine Pyrophosphate
1256. Thiamine Triphosphate
1257. Thiamylal
1258. Thienylcyclohexylpiperidine
1259. Thiobarbituric Acid
1260. Thiochrome
1261. Thiochrome Monophosphate
1262. Thiochrome Pyrophosphate
1263. Thiochrome Triphosphate
1264. Thiocyanate
1265. Thionin
1266. Thiosalicylic Acid
1267. Thiosulfate
1268. Thiourea
1269. Threonine
1270. Thymidine
1271. Thymidine 5'-diphosphate (TDP)
1272. Thymidine 5'-monophosphate (TMP)
1273. Thymidine 5'-triphosphate (TTP)
1274. Thymine
1275. Thymol (2-Isopropyl-5-Methylphenol)
1276. Thymyl Methyl Ether
1277. Thyrotropin-releasing Hormone
1278. Thyrotropin-releasing Hormone Metabolites
1279. Thyroxine (T4)
1280. Tinopal 5BM-GX
1281. Tinopal CBS-X
1282. Tin Protoporphyrin
1283. Tocopherol
1284. <i>α</i> -Tocopherol
1285. <i>d,l</i> - <i>α</i> -Tocopherol
1286. <i>g</i> -Tocopherol
1287. Toluene
1288. <i>o</i> -Toluidine
1289. Transferrin
1290. Trehalose
1291. Triallylsulfonium
1292. Tributylamine
1293. Tricaprin
1294. 1,2,3-Trichlorobenzene
1295. 1,2,4-Trichlorobenzene
1296. 1,2,5-Trichlorobenzene
1297. 2,3,4-Trichlorophenol

1298.	2,3,5-Trichlorophenol
1299.	2,3,6-Trichlorophenol
1300.	2,4,5-Trichlorophenol
1301.	2,4,6-Trichlorophenol
1302.	3,4,5-Trichlorophenol
1303.	2,4,5-Trichlorophenoxyacetic Acid
1304.	(+)-2-(2,4,5-Trichlorophenoxy) propionic Acid
1305.	Triclocarban
1306.	Tridecane
1307.	Triethanolamine
1308.	Triethylamine
1309.	Trifluoroacetic Acid
1310.	1,2,3-Trihydroxybenzene
1311.	Trimethylamine
1312.	Trimethylarsine Oxide
1313.	Trimethylselenonium Iodide
1314.	Trimethylselenonium Ion
1315.	Trimethylsulfonium Ion
1316.	Triphenyltin Hydroxide
1317.	Tripolyphosphate
1318.	Tris(Acetylacetonato)Aluminum (III)
1319.	Tris(Acetylacetonato)Chromium (III)
1320.	Tris(Acetylacetonato)Iron (III)
1321.	Tris Buffer
1322.	Tris(Hexafluoroacetylacetonato) Chromium (III)

1323.	Tris(1-Phenyl-1,3-Butane-Dionato) Chromium (III)
1324.	Tromethamine
1325.	Trypsin
1326.	Tryptic Digest of Cytochrome c
1327.	Tryptophan
1328.	Tyrosine
1329.	Tween 80

U, V, W

1330.	Undecane
1331.	1-Undecanol
1332.	2-Undecanone
1333.	Uracil
1334.	Urethane
1335.	Uric Acid
1336.	Uridine
1337.	Uridine 5'-Monophosphate (UMP)
1338.	Ursodeoxycholate, Glycine and Taurine-conjugated
1339.	Ursodeoxycholic Acid
1340.	Valerate
1341.	Valeric Acid
1342.	a-Keto, b-Methyl, n-Valeric Acid
1343.	Valerophenone
1344.	Valine
1345.	Vanillic Acid

1346.	Vanillin
1347.	Verdamicin
1348.	Vitamin A
1349.	Vitamin B12
1350.	Vitamin D2
1351.	Vitamin E
1352.	Vitamin K1
1353.	Warfarin
1354.	Water

X, Y, Z

1355.	Xanthine
1356.	Xanthosine
1357.	2,3-Xylenol
1358.	2,4-Xylenol
1359.	2,5-Xylenol
1360.	2,6-Xylenol
1361.	3,4-Xylenol
1362.	3,5-Xylenol
1363.	Xylitol
1364.	Xylose
1365.	Yohimbine
1366.	Zinc
1367.	Zinc Deuteroporphyrin
1368.	Zinc Phenanthroline
1369.	Zinc Protoporphyrin

Can't find your compound inside? Call us at 800-648-5950 ext. 290 or e-mail: hplc@hamiltoncompany.com. We'll tell you if your compound can be separated before you buy a column, if we can't separate your compounds you haven't spent a dime.

HAMILTON
THE MEASURE OF EXCELLENCE™

Hamilton Company
P.O. Box 10030
Reno, Nevada 89520-0012 U.S.A.
Toll-Free 800-648-5950
Fax +1-775-856-7259
Telephone +1-775-858-3000
e-mail: hplc@hamiltoncompany.com
<http://www.hamiltoncompany.com>

