

Crimper and Decapper



Easy Grip Manual Crimper & Decapper
Electronic Crimper & Decapper
High Power Electronic Crimper & Decapper
Classic Manual Crimper & Decapper

BGB Easy Grip Manual Crimper and Decapper

Our Easy Grip Manual Crimper and Decapper for all 8, 11, 13 and 20 mm crimp caps. The curved handle features comfortable use compared to metal grip designs. The bottom-pull handle allows for a steadier hold than top-push handles.

An easily-viewed adjustment knob is located at the head of the crimper and decapper. The clearly shown + and - symbols along with the directional arrows, simplify adjusting the intensity of the crimp or decap needed. The knob is also designed to indicate when the crimp setting has been reached, when it fully contacts the body of the crimper/decapper.



All Easy Grip Manual Crimper and Decapper are equipped with an adjustment knob.



Part No.	Description	Qty.
295000	8mm Crimper	ea.
295001	11 mm Crimper	ea.
295002	13 mm Crimper	ea.
295003	13 mm Crimper for Flip Off and Tear Off Caps	ea.
295004	20 mm Crimper	ea.
295005	20 mm Crimper for Flip Off and Tear Off Caps	ea.
295006	11 mm Decapper	ea.
295007	13 mm Decapper	ea.
295008	20 mm Decapper	ea.



Easy to use.



BGB Classic Manual Crimper and Decapper

- Chemically resistant surface finish especially designed for the application in the labs.
- Hardened crimping jaws made of special alloy that guarantees long life
- Adjustable in crimping pressure limitation by a screw in the handle
- Additionally 11 mm, 13 mm and 20 mm crimpers are adjustable in crimping height by screwing up or down the pressure block in the crimping head with a hexagon key.



Part No.	Description	Qty.
295020	8mm Crimper	ea.
295021	11 mm Adjustable Crimper	ea.
295022	13 mm Adjustable Crimper	ea.
295023	13 mm Adjustable Crimper for Flip Off and Tear Off Caps	ea.
295024	20 mm Adjustable Crimper	ea.
295025	20 mm Adjustable Crimper or Flip Off and Tear Off Caps	ea.
295026	8mm Decapper	ea.
295027	11 mm Decapper	ea.
295028	13 mm Decapper	ea.
295029	20 mm Decapper	ea.

BGB Electronic Crimper and Decapper Model 6 New

Our new Electronic Crimpers and Decappers with enhanced OLED screen displays fit on standard laboratory vials and seals with aluminum caps. The electronically adjustable cycle is accurate and reproducible and can be adjusted on the display of the unit. The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools.

- Powerful motor with fast cycle
- Longer lasting battery with Lithium Ion technology
- Advanced charging system including warning light when recharging is necessary
- Can be operated while charging
- Controls on top of unit with enhanced OLED Screen



OLED screen display



Easy to use.



Part No.	Description	Qty.
297040	8 mm Crimper Model 6	ea.
297041	11 mm Crimper Model 6	ea.
297042	13 mm Crimper Model 6	ea.
297043	13 mm Crimper Model 6 for Flip Off and Tear Off Caps	ea.
297044	20 mm Crimper Model 6	ea.
297045	20 mm Crimper Model 6 for Flip Off and Tear Off Caps	ea.
297046	11 mm Decapper Model 6	ea.
297047	13 mm Decapper Model 6	ea.
297048	20 mm Decapper Model 6	ea.

Note: Not suitable for all-steel or magnetic ring crimp caps, please use High Power Electronic Crimping Tool (see page 4) for this type of caps.

BGB Crimping Tool Stand

Stand for BGB Electronic and Easy Grip Manual Crimpers and Decappers.



Part No.	Description	Qty.
295080	Crimping Tool Stand	ea.



Examples of use.

BGB High Power Electronic Crimping Tool HPS Model 6



The High Power Electronic Crimping Tool with new OLED screen display is faster and more powerful than the electronic crimping tool. The HPS crimping tool is designed to store separate programs for different caps and seals. It is powered by a universal power supply and uses interchangeable jaw sets. Jaw sets are strong enough to crimp all steel and magnetic caps.

- High-Power version of electronic crimping tool with fixed power supply and cord (no battery).
- Uses interchangeable jaw sets.
- Stores separate programs for different caps and seals.
- Strong enough for all steel and magnetic caps.
- Optional base with mounting kit.



OLED screen display

Part No.	Description	Qty.
297060	High Power Electr. Crimping Tool Model 6	ea.



High Power Electronic Crimping Tool is not included.

Part No.	Description	Qty.
295061	Base with Mounting Kit	ea.

Jaw Sets for BGB High Power Electronic Crimping Tool HPS

Dimension Tool	8 mm Crimper	11 mm Crimper	13 mm Crimper	13 mm Crimper for Flip Off and Tear Off Caps	20 mm Crimper	20 mm Crimper for Flip Off and Tear Off Caps
Part No.	295062	295063	295064	295065	295066	295067
Qty.	ea.	ea.	ea.	ea.	ea.	ea.

Dimension Tool	11 mm Decapper	13 mm Decapper	20 mm Decapper
Part No.	295068	295069	295070
Qty.	ea.	ea.	ea.



Jaw Sets are easy to change. Simply insert the jaw set into the bushing at the bottom of the tool. Push up and twist until the set locks into position. To remove, push the outside button and rotate.

Prices may change without notice and do not include VAT.