

Product Overview

- Adjustable height from 57 cm to 87 cm
- 2 liftable columns
- Trespa™ chemical resistant worksurface W90 x D75
- Weight capacity: 250 kg
- Worksurface can easily be lifted up or down by commutating an electrical switch
- Worldwide electrical compatibility
- Recommended use with liquid chromatograph
- HPLC Ionbench is on caster wheels and can easily be moved
- Convenient access to the top of your HPLC or UHPLC
- Safe access to solvents

Liftable
motorized
worksurface



Control :



Standard switch

- *Up and Down*



Optional position display with 3 memories (ref : BCHPOS)

- *Up and Down*
- *Memory function (3 in total)*
- *Display for height indication*

Photo example :



Detailed Technical specifications :

- Compatible with 230V/50HZ or 110V/60HZ
- Electronic components are all compliant with RoHS directives (Restriction of Hazardous Substances Directive).

External Dimensions	BCHLC9075
Depth	75 cm / 29,5 inch
Width	90 cm / 35,4 inch
Height	From 57 to 87 cm / From 22,4 to 34,3 inch
Weight :	35 kg / 77 lb

Summary

Description	Part-number
Ionbench-LC BCHLC9075 / 220 Volt	BCHLC9075V220
Ionbench-LC BCHLC9075 / 110 Volt	BCHLC9075V110
Optional position display with 3 memories	BCHPOS