


Inlet Supplies for PerkinElmer GCs With CAP Inlet



Septa (11 mm) for PerkinElmer GCs

Thermolite Septa	
Usable to 340 °C inlet temp.	
qty.	cat.#
50-pk.	27142
100-pk.	27143

Premium Non-Stick BTO Septa	
Usable to 400 °C inlet temp.	
qty.	cat.#
50-pk.	27090
100-pk.	27091

Silicone O-Rings

Description	Max. temp.	Similar to PE part #	qty.	cat.#
Silicone O-Rings for PerkinElmer Auto SYS XL or Clarus w/CAP Injector	250 °C	N6101374	10-pk.	20262

Sky® Inlet Liner

ID x OD x L	qty.	cat.#
Precision, Sky Technology, Borosilicate Glass with Quartz Wool	ea.	23450.1
4.0 mm x 6.2 mm x 92.1 mm	5-pk.	23450.5

Capillary Ferrules for 1/16-Inch Compression-Type Fittings

Ferrule ID	Fits Column ID	qty.	Vespel/Graphite (60/40)
0.4 mm	0.10/0.15/0.18/0.25/0.28 mm	10-pk.	20211
0.5 mm	0.32 mm	10-pk.	20212
0.8 mm	0.45/0.53 mm	10-pk.	20213

RESTEK
Pure Chromatography