




## Fast separations on short CHIRAL-CBH columns

Due to the extremely high enantioselectivity it is easy to obtain very fast separations on short **CHIRAL-CBH** columns.

Below are some examples on compounds that have been resolved with retention times for the second eluted peak under 2.2 min.



**CHIRAL-CBH** columns are available in the following analytical dimensions:

50x2.0 mm    50x3.0 mm    50x4.0 mm  
100x2.0 mm    100x3.0 mm    100x4.0 mm  
150x2.0 mm    150x3.0 mm    150x4.0 mm

Guard cartridges matching the above columns (holder necessary):

10x2 mm    10x3 mm    10x4 mm

Semipreparative columns with ID 10 mm are available in lengths of 100 and 150 mm.

We also supply other dimensions on request, for example capillaries down to ID 50 um.

Many more applications on [www.chromtech.co.uk](http://www.chromtech.co.uk)